

Product Datasheet

ACD Antibody (1C11-1A7) H00065057-M01

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00065057-M01

Updated 5/13/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00065057-M01



H00065057-M01


ACD Antibody (1C11-1A7)

Product Information	
Unit Size	0.1 mg
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	1C11-1A7
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Ascites
Buffer	In 1x PBS, pH 7.4
Product Description	
Host	Mouse
Gene ID	65057
Gene Symbol	ACD
Species	Human
Specificity/Sensitivity	ACD - adrenocortical dysplasia homolog (mouse) (1C11-1A7)
Immunogen	ACD (AAH16904, 1 a.a. ~ 544 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa. MPGRCQSDAAMRVNGPASRAPAGWTSGSLHTGPRAGRPRQAARGVRGRGL LLRPRPAKELPLPRKGGAWAPAGNPGPLHPLGVAVGMAGSGRLVLRPWIRELI LGSETPSSPRAGQLLEVLQDAEAAVAGPSHAPDTSVDGATLLVSDGTHSVRCL VTREALDTSDWEEKEFGFRGTEGRLLLLQDCGVHVQVAEGGAPAEFYLVQVDR FSLLPTEQPRLRVPGCNQDLVDVQKKLYDCLEEHLESTSSNAGLSLSQLLDEM REDQEHQALVCLAESCLTLEGPCTAPPVTHWAASRCKATGEAVYTPSSML CISENDQLILSSLGPCQRTQGPELPPDPALQDLSLTIASPPSSPSSSGTPALP GHMSSEESGTSISLLPALSLAAPDPGQRSSSQPSAICSAPATLTTPRSPHASRT PSSPLQSCTPSLSPRSHVPSPHQALVTRPQKPSLEFKEFVGLPCKNRPPFPRT GATRGAEPCSVWEPPKRRHRDGSFQYEEPPCTSLCARVQAARLPPQLMA WALHFLMDAQPGSEPTPM
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, ELISA, Immunoprecipitation, Sandwich ELISA
Recommended Dilutions	Western Blot 1:100 - 1:2000, ELISA 1:100 - 1:2000, Immunoprecipitation 1:10 - 1:500, Sandwich ELISA




Images


Western Blot: ACD Antibody (1C11-1A7) [H00065057-M01] - Analysis of ACD expression in transfected 293T cell line by ACD monoclonal antibody (M01), clone 1C11-1A7. Lane 1: ACD transfected lysate (Predicted MW: 57.7 KDa). Lane 2: Non-transfected lysate.



Immunoprecipitation: ACD Antibody (1C11-1A7) [H00065057-M01] - Analysis of ACD transfected lysate using anti-ACD monoclonal antibody and Protein A Magnetic Bead, and immunoblotted with ACD MaxPab rabbit polyclonal antibody.



Sandwich ELISA: ACD Antibody (1C11-1A7) [H00065057-M01] - Detection limit for recombinant GST tagged ACD is 1 ng/ml as a capture antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00065057-M01

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

