

# Product Datasheet

## DSCC1 293T Cell Transient Overexpression Lysate H00079075-T01

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00079075-T01](http://www.novusbio.com/H00079075-T01)

Updated 6/27/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00079075-T01](http://www.novusbio.com/reviews/destination/H00079075-T01)



**H00079075-T01****DSCC1 293T Cell Transient Overexpression Lysate**

| <b>Product Information</b> |                                                                                                           |
|----------------------------|-----------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>           | 0.1 ml                                                                                                    |
| <b>Concentration</b>       | Please contact technical services for concentration.                                                      |
| <b>Storage</b>             | Store at -80C. Avoid freeze-thaw cycles.                                                                  |
| <b>Buffer</b>              | 1X Sample Buffer (50 mM Tris-HCl, 2% SDS, 10% glycerol, 300 mM 2-mercaptoethanol, 0.01% Bromophenol blue) |

| <b>Product Description</b>     |                                                                                                                                                       |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>             | Quality control test: Transient overexpression cell lysate was tested with Anti-DCC1 antibody by Western Blots. <b>Plasmid:</b> pCMV-DCC1 full-length |
| <b>Gene ID</b>                 | 79075                                                                                                                                                 |
| <b>Gene Symbol</b>             | DSCC1                                                                                                                                                 |
| <b>Species</b>                 | Human                                                                                                                                                 |
| <b>Specificity/Sensitivity</b> | DCC1 293T Cell Transient Overexpression Lysate(Denatured)                                                                                             |
| <b>Notes</b>                   | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                    |
| <b>Lysate Type</b>             | 293T Cell Transient Overexpression                                                                                                                    |
| <b>Lysate Protein State</b>    | Denatured                                                                                                                                             |

| <b>Product Application Details</b> |              |
|------------------------------------|--------------|
| <b>Applications</b>                | Western Blot |
| <b>Recommended Dilutions</b>       | Western Blot |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to H00079075-T01**

---

|                   |                                                    |
|-------------------|----------------------------------------------------|
| H00079075-P01-2ug | Recombinant Human DSCC1 GST (N-Term) Protein       |
| AF2535            | MMR/CD206/Mannose Receptor Antibody [Unconjugated] |
| H00079075-B01P    | DSCC1 Antibody                                     |
| NB500-106         | PCNA Antibody (PC10)                               |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00079075-T01](http://www.novusbio.com/reviews/submit/H00079075-T01)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

