

Product Datasheet

FLJ21908 Antibody H00079657-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00079657-B01P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00079657-B01P



H00079657-B01P

FLJ21908 Antibody

Product Information

Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)

Product Description

Host	Mouse
Gene ID	79657
Gene Symbol	RPAP3
Species	Human
Specificity/Sensitivity	RPAP3 - RNA polymerase II associated protein 3,
Immunogen	RPAP3 (AAH56415.1, 1 a.a. - 631 a.a.) full-length human protein. MTSANKAIELQLQVKQNAEELQDFMRDLENWEKDIKQKDMELRRQNGVPEEN LPPIRNGNFRKKKKGKAKESSKKTREENTKNRIKSVDYEAWAKLDVDRILDELD KDDSTHESLSQESESEEDGIHVDSQKALVLKEKGKNGYFKQKGYDEAIDCYTKG MDADPYNPVLPTNRASAYFRLKKFAVAESDCNLAVALNRSYTKAYSRRGAARF ALQKLEEAKKDYERVLELEPNNFEATNELRKISQALASKENSYPKEADIVIKSTE GERKQIEAQNKQQAISEKDRGNGFFKEGKYERAIECYTRGIAADGANALLPA NRAMAYLKIQKYEEAEKDCTQAILLDGSYSKAFARRGTARTFLGKLNEAKQDFE TVLLLEPGNKQAVTELSKIKKKPLKKVIIIEETGNLIQTIDVPDSTTAAAPENNPINL ANVIAATGTTSKKNSSQDDLFPSTDPRAKVLKIEEVSDTSSLQPQASLKQDVC QSYSEKMPIEIEQKPAQFATTVLPIIPANSFQLESDFRQLKSSPDMLYQYLKQIE PSLYPKLFQKNLDPDVFNQIVKILHDFYIEKEKPLLIFEILQRLSELKRFDMAVMF MSETEKKIARALFNHIDKSGLDSSVEELKKRYGG
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.


Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1:500
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.



Images

Western Blot: FLJ21908 Antibody [H00079657-B01P] - Analysis of RPAP3 expression in transfected 293T cell line by RPAP3 polyclonal antibody. Lane 1: FLJ21908 transfected lysate(69.41 KDa). Lane 2: Non-transfected lysate.



Publications

Izumi N, Yamashita A, Hirano H, Ohno S. Heat shock protein 90 regulates phosphatidylinositol 3-kinase-related protein kinase family proteins together with the RUVBL1/2 and Tel2-containing co-factor complex. Cancer Sci. 2012-01-01 [PMID: 21951644] (Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00079657-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

