

# Product Datasheet

## PRKRIP1 Antibody (3B11) - Azide and BSA Free H00079706-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00079706-M02](http://www.novusbio.com/H00079706-M02)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00079706-M02](http://www.novusbio.com/reviews/destination/H00079706-M02)




**H00079706-M02****PRKRIP1 Antibody (3B11) - Azide and BSA Free**

| <b>Product Information</b>         |                                                                                                                                                                                                                                                                                                                                |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>                   | 0.1 mg                                                                                                                                                                                                                                                                                                                         |
| <b>Concentration</b>               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                               |
| <b>Storage</b>                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                                                                                                                                                                                                                   |
| <b>Clonality</b>                   | Monoclonal                                                                                                                                                                                                                                                                                                                     |
| <b>Clone</b>                       | 3B11                                                                                                                                                                                                                                                                                                                           |
| <b>Preservative</b>                | No Preservative                                                                                                                                                                                                                                                                                                                |
| <b>Isotype</b>                     | IgG1 Kappa                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>                      | IgG purified                                                                                                                                                                                                                                                                                                                   |
| <b>Buffer</b>                      | In 1x PBS, pH 7.4                                                                                                                                                                                                                                                                                                              |
| <b>Product Description</b>         |                                                                                                                                                                                                                                                                                                                                |
| <b>Description</b>                 | Novus Biologicals Mouse PRKRIP1 Antibody (3B11) - Azide and BSA Free (H00079706-M02) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                           |
| <b>Host</b>                        | Mouse                                                                                                                                                                                                                                                                                                                          |
| <b>Gene ID</b>                     | 79706                                                                                                                                                                                                                                                                                                                          |
| <b>Gene Symbol</b>                 | PRKRIP1                                                                                                                                                                                                                                                                                                                        |
| <b>Species</b>                     | Human                                                                                                                                                                                                                                                                                                                          |
| <b>Reactivity Notes</b>            | Human. Other species not tested.                                                                                                                                                                                                                                                                                               |
| <b>Specificity/Sensitivity</b>     | PRKRIP1 - PRKR interacting protein 1 (IL11 inducible)                                                                                                                                                                                                                                                                          |
| <b>Immunogen</b>                   | PRKRIP1 (AAH14298, 1 a.a. ~ 184 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.<br>MASPAASSVRPPRPKKEPQTLVIPKNAAEEQKLKLERLMKNPDKAVPIPEKMSE<br>WAPRPPPEFVRDVMGSSAGAGSGEFHVYRHLRRREYQRQDYMDAMAEKQK<br>LYAEFQKRLEKNKIAAEEQTAKRRKKRQKLKEKKLLAKKMKLEQKKQEGPGQP<br>KEQGSSSSAEASGTEEEEEEVPSFTMGR |
| <b>Notes</b>                       | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                             |
| <b>Product Application Details</b> |                                                                                                                                                                                                                                                                                                                                |
| <b>Applications</b>                | Western Blot, ELISA                                                                                                                                                                                                                                                                                                            |
| <b>Recommended Dilutions</b>       | Western Blot 1:500, ELISA                                                                                                                                                                                                                                                                                                      |
| <b>Application Notes</b>           | Antibody Reactive against Recombinant Protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control.                                                                                                                                                                                             |



## Images

Western Blot: PRKRIP1 Antibody (3B11) [H00079706-M02] - Detection against Immunogen (45.98 KDa).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00079706-M02**

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-43319-0.5mg   | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)           |
| H00079706-P01-10ug | Recombinant Human PRKRIP1 GST (N-Term) Protein          |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00079706-M02](http://www.novusbio.com/reviews/submit/H00079706-M02)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

