

# Product Datasheet

## MOGAT2 Antibody - Azide and BSA Free H00080168-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00080168-B01P](http://www.novusbio.com/H00080168-B01P)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00080168-B01P](http://www.novusbio.com/reviews/destination/H00080168-B01P)



**H00080168-B01P**

MOGAT2 Antibody - Azide and BSA Free

**Product Information**

|                      |                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.05 mg                                                                                                          |
| <b>Concentration</b> | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| <b>Clonality</b>     | Polyclonal                                                                                                       |
| <b>Preservative</b>  | No Preservative                                                                                                  |
| <b>Isotype</b>       | IgG                                                                                                              |
| <b>Purity</b>        | Protein G purified                                                                                               |
| <b>Buffer</b>        | PBS (pH 7.4)                                                                                                     |

**Product Description**


|                    |                                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b> | Novus Biologicals Mouse MOGAT2 Antibody - Azide and BSA Free (H00080168-B01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                                                                                                   |
| <b>Host</b>        | Mouse                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Gene ID</b>     | 80168                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Gene Symbol</b> | MOGAT2                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Species</b>     | Human                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Immunogen</b>   | MOGAT2 (AAI03879.1, 1 a.a. - 284 a.a.) full-length human protein.<br>MVEFAPLFMPWERRLQTLAVLQFVFSFLALAEICTVGFIALLFTRFWLLTVLYAA<br>WWYLDNRDKPRQGGRHIQAIRCWTIWKYMKDYPISLVKTAELDPSRNYIAGFH<br>PHGVLAVGAFANLCTESTGFSSIFPGIRPHLMMLTLWFRAPFFRDYIMSAGLVT<br>SEKESAAHILNRKGGGNLLGIIVGGAQEALDARPGSFTLLLRNRKGFVRLALTH<br>GYQASGKSTLGSVGNWQGFYFGGKMAETNADSILVEIFSPFTIKIIFWCLMPKY<br>LEKFPQRRLSDLRN |
| <b>Notes</b>       | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                                                                                    |

**Product Application Details**

|                              |                                                      |
|------------------------------|------------------------------------------------------|
| <b>Applications</b>          | Western Blot                                         |
| <b>Recommended Dilutions</b> | Western Blot                                         |
| <b>Application Notes</b>     | This antibody is useful for Western Blot, Functional |

**Images**

Western Blot: MOGAT2 Antibody [H00080168-B01P] - Analysis of MOGAT2 expression in transfected 293T cell line by MOGAT2 polyclonal antibody. Lane 1: MOGAT2 transfected lysate(31.24 KDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

[nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

[canada.inquires@bio-techne.com](mailto:canada.inquires@bio-techne.com)

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

[info.EMEA@bio-techne.com](mailto:info.EMEA@bio-techne.com)

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)

Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)

Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)

General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to H00080168-B01P**

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97019-5mg     | Mouse IgG Isotype Control                               |
| H00080168-P01-10ug | Recombinant Human MOGAT2 GST (N-Term) Protein           |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00080168-B01P](http://www.novusbio.com/reviews/submit/H00080168-B01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

