

Product Datasheet

ILKAP/PP2C delta Antibody

H00080895-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00080895-B01P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00080895-B01P



H00080895-B01P

ILKAP/PP2C delta Antibody

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Mouse
Gene ID	80895
Gene Symbol	ILKAP
Species	Human
Specificity/Sensitivity	ILKAP - integrin-linked kinase-associated serine/threonine phosphatase 2C,
Immunogen	ILKAP (NP_110395.1, 1 a.a. - 392 a.a.) full-length human protein. MDLFGDLPEPERSPRPAAGKEAQKGPLLFDLPPASSTDGSGGGPLLFDLPP ASSGDSGLATSISQMVKTEGKGAKRKTSEEEKNGSEELVEKKVCKASSVIFG LKG YVAERKGEREEMQDAHVILNDITEECRPPSSLITRVSYFAVFDGHGGIRAS KFAAQNLHQNLIRKFPKGDVISVEKTVKRCLLDTFKHTDEEFLKQASSQKPAWK DGSTATCVLAVDNILYIANLGDSRAILCRYNEESQKHAALSLSKEHNPTQYEER MRIQKAGGNVRDGRVLGVLEVSRSIGDGQYKRCGVTSVPDIRRCQLTPNDRFI LLACDGLFKVFTPEEAVNFILSCLEDEKIQTREGKSAADARYEAACNRLANKAV QRGSADNVTVMVVRIGH
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/Immunofluorescence
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.



Images

Western Blot: ILKAP/PP2C delta Antibody [H00080895-B01P] - Analysis of ILKAP expression in transfected 293T cell line by ILKAP polyclonal antibody. Lane 1: ILKAP transfected lysate(43.12 KDa). Lane 2: Non-transfected lysate.



Publications

Zhou W, Cao H, Yang X et al. Characterization of Nuclear Localization Signal in the N Terminus of ILKAP and its Essential Role in the Down-Regulation of RSK2 Signaling. J Biol Chem. 2013-01-17 [PMID: 23329845]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00080895-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

