

# Product Datasheet

## STK22C Antibody - Azide and BSA Free H00081629-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00081629-D01P](http://www.novusbio.com/H00081629-D01P)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00081629-D01P](http://www.novusbio.com/reviews/destination/H00081629-D01P)



**H00081629-D01P**

STK22C Antibody - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Rabbit STK22C Antibody - Azide and BSA Free (H00081629-D01P) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-STK22C Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	81629
Gene Symbol	TSSK3
Species	Human, Mouse
Specificity/Sensitivity	TSSK3 - testis-specific serine kinase 3,
Immunogen	TSSK3 (NP_443073.1, 1 a.a. - 268 a.a.) full-length human protein. MEDFLLSNGYQLGKTIGEGTYSKVKEAFSKKHQRKVAIKVIDKMGGPEEFIQRFLPRELQIVRTL DHKNIIQVYEMLESADGKICLVMELAEAGGDVFDCVLNNGGPLPESRAKALFRQMVEAIRYCHGCGVAHRDLKCENALLQG FNKL TDFGFAKVLPKSHRELSQTFCGSTAYAAPEVLQGIPHDSKKGDVWSMGVVLYVMLCASLPFDDTDIPKMLWQQQKGVSFPTHLSISADCQDLLKRLLEPDMILRPSIEEVSWHPWLAS T
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control. Immunocytochemistry/Immunofluorescence and Immunohistochemistry were reported in scientific literature.




## Images

Western Blot: STK22C Antibody [H00081629-D01P] - Analysis of TSSK3 expression in transfected 293T cell line by TSSK3 polyclonal antibody. Lane 1: TSSK3 transfected lysate(30.10 KDa). Lane 2: Non-transfected lysate.



Western Blot: STK22C Antibody [H00081629-D01P] - TSSK3 MaxPab rabbit polyclonal antibody. Western Blot analysis of TSSK3 expression in mouse liver.



## Publications

Li Y, Sosnik J, Brassard L et al. Expression and Localization of Five Members of the Testis-Specific Serine Kinase (Tssk) Family in Mouse and Human Sperm and Testis. Mol Hum Reprod 2010-08-20 [PMID: 20729278] (WB, ICC/IF, IF/IHC, Mouse)

**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to H00081629-D01P**

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00081629-P01-10ug	Recombinant Human STK22C GST (N-Term) Protein

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00081629-D01P](http://www.novusbio.com/reviews/submit/H00081629-D01P)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

