

# Product Datasheet

## ARPC5L Overexpression Lysate (Denatured) H00081873-T01

Unit Size: 0.1 ml

Store at -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00081873-T01](http://www.novusbio.com/H00081873-T01)

Updated 8/18/2016 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00081873-T01](http://www.novusbio.com/reviews/destination/H00081873-T01)




**H00081873-T01****ARPC5L Overexpression Lysate (Denatured)**


| <b>Product Information</b>         |  |
|------------------------------------|--|
| <b>Unit Size</b>                   | 0.1 ml   |
| <b>Concentration</b>               | The exact concentration of the protein of interest cannot be determined for overexpression lysates. Please contact technical support for more information. |
| <b>Storage</b>                     | Store at -80C. Avoid freeze-thaw cycles.   |
| <b>Buffer</b>                      | 50mM Tris-HCl, 2%SDS, 10% glycerol, 2% 2-mercaptoethanol, pH 6.8.  |
| <b>Product Description</b>         |  |
| <b>Gene ID</b>                     | 81873  |
| <b>Gene Symbol</b>                 | ARPC5L   |
| <b>Species</b>                     | Human  |
| <b>Reactivity Notes</b>            | Human  |
| <b>Specificity/Sensitivity</b>     | actin related protein 2/3 complex, subunit 5-like  |
| <b>Notes</b>                       | This product is produced by and distributed for Abnova, a company based in Taiwan.   |
| <b>Lysate Type</b>                 | Overexpression   |
| <b>Lysate Protein State</b>        | Denatured  |
| <b>Product Application Details</b> |  |
| <b>Applications</b>                | Western Blot   |
| <b>Recommended Dilutions</b>       | Western Blot   |
| <b>Application Notes</b>           | This transient overexpression cell lysate is useful for western blot and was tested with Anti-ARPC5L antibody H00081873-B01 by Western Blots.              |

**Images**

ARPC5L transfected lysate.



Lane 1: ARPC5L transfected lysate (16.94 KDa) Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

8100 Southpark Way, A-8  
Littleton, CO 80120  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Novus Biologicals Canada**

461 North Service Road West, Unit B37  
Oakville, ON L6M 2V5  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada@novusbio.com

### **Novus Biologicals Europe**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to H00081873-T01**

---

|                    |                            |
|--------------------|----------------------------|
| H00081873-P01-10ug | ARPC5L Recombinant Protein |
| NBP2-15472         | ARPC5L Antibody            |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Lysates are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00081873-T01](http://www.novusbio.com/reviews/submit/H00081873-T01)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

