

Product Datasheet

FERMT3/URP2 Antibody H00083706-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00083706-B01P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00083706-B01P



H00083706-B01P**FERMT3/URP2 Antibody**

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	IgG purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Mouse
Gene ID	83706
Gene Symbol	FERMT3
Species	Human
Specificity/Sensitivity	URP2 - UNC-112 related protein 2,
Immunogen	FERMT3 (NP_113659.3, 1 a.a. - 663 a.a.) full-length human protein. MAGMKTASGDYIDSSWELRVFVGEEDPEAESVTLRVTGESHIGGVLLKIVEQIN RKQDWSHDHAIWWEQKRQWLLQTHWTLDKYGILADARLFFGPQHRPVILRLPN RRALRLRASFSQPLFQAVAAICRLLSIRHPEELSLLRAPEKKEKKKKEKEPEEEL YDLSKVVLGGVAPALFRGMPAHFSDSAQTEACYHMLSRPQPPPDPLLLQRLP RPSSLSDKTQLHSRWLDSSRCLMQQGKAGDALWLRFKYYSFFDLDPKTDPV RLTQLYEQARWDLLEEDCTEEEMMVFAALQYHINKLSQSGEVGEPAAGTDPG LDDLVALSNLEVKLEGSAPTDVLDLSTTIPELKDHLRIFRPRKLTCLKGYRQHWV VFKETTLSYYKSQDEAPGDPQQNLKKGCEVVPDVNVSGQKFCIKLLVPSPEG MSEIYLRCQDEQQYARWMAGCRLASKGRMTADSSYTSEVQAILAFLSLQRTG SGGPGNHPHGPDASAEGLNPYGLVAPRFQRKFKAKQLTPRILEAHQNVAQLSL AEAQLRFIQAWQSLPDFGISYVMVRFKGSRKDEILGIANNRLRIDLAVGDVVKT WRFSNMRQWNVNWDIRQVAIEFDEHINVAFCVSVASCRIVHEYIGGYIFLSTRE RARGEELDEDLFLQLTGGHEAF
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/Immunofluorescence
Application Notes	Antibody reactive against transfected lysate for Western Blot. Has also been used for immunofluorescence and ELISA.




Images

Western Blot: FERMT3/URP2 Antibody [H00083706-B01P] - Analysis of FERMT3 expression in transfected 293T cell line by FERMT3 polyclonal antibody. Lane 1: URP2 transfected lysate(72.93 KDa). Lane 2: Non-transfected lysate.



Immunocytochemistry/Immunofluorescence: FERMT3/URP2 Antibody [H00083706-B01P] - Analysis of purified antibody to URP2 on HeLa cell. (antibody concentration 10 ug/ml)



Publications

Feng C, Li YF, Yau YH et al. Kindlin-3 mediates integrin alphaLb2 outside-in signaling and it interacts with the scaffold protein receptor for activated-C kinase 1 (RACK1). J Biol Chem. 2012-02-10 [PMID: 22334666]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00083706-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

