

Product Datasheet

SNX25 Antibody - Azide and BSA Free H00083891-D01P

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00083891-D01P

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00083891-D01P



H00083891-D01P

SNX25 Antibody - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)


Product Description	
Description	Novus Biologicals Rabbit SNX25 Antibody - Azide and BSA Free (H00083891-D01P) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	83891
Gene Symbol	SNX25
Species	Human
Specificity/Sensitivity	SNX25 - sorting nexin 25,
Immunogen	SNX25 (ENSP00000296779, 1 a.a. - 483 a.a.) full-length human protein. MYMERMDKRALISFWESVEHLKNANKNEIPQLVGEIYQNFFVESKEISVEKSLY KEIQQCLVGNKGIEVFYKIQEDVYETLKDRYYPSFIVSDLYEKLLIKEEEKHASQ MISNKDEMGRDEAGEEAVDDGTNQINEQASFAVNKLRELNEKLEYKRQALNS IQNAPKPKDKIVSKLKDEIILIEKERTDLQLHMARTDWWCENLGMWKASITSGEV TEENGEQLPCYFVMVSLQEVGGVETKNWTVPRRLSEFQNLHRKLSECVPSLK KVQLPSLSKLPFKSIDQKFMEKSKNQLNKFLQEETEEDSDLSYDGDVDGRKD ALAEPCFMLIGEIFELRGMFKWVRRTLIALVQVTFGRTINKQIRD TVSWIFSEQM LVYYINIFRDAFWPNGKLAPPTTIRSKEQSQETKQRAQQKLLLENIPDMLQSLVG QQNARHGIKIFNALQETRANKHLLYALMELLIELCPELRVHLDQLKAGQV
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1:500
Application Notes	Antibody reactive against Recombinant Protein with GST tag on ELISA and Western Blot and also on transfected lysate in western blot. GST tag alone is used as a negative control.



Images

Western Blot: SNX25 Antibody [H00083891-D01P] - Analysis of SNX25 expression in transfected 293T cell line by SNX25 polyclonal antibody. Lane 1: SNX25 transfected lysate(56.30 KDa). Lane 2: Non-transfected lysate.



**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to H00083891-D01P

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00083891-P01-10ug	Recombinant Human SNX25 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00083891-D01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

