

Product Datasheet

Ln timer Antibody **H00084708-B01P**

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00084708-B01P

Updated 1/12/2021 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00084708-B01P



H00084708-B01P


Lnx1 Antibody

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	84708
Gene Symbol	LNK1
Species	Human
Immunogen	LNK1 (AAH34737.1, 1 a.a. - 632 a.a.) full-length human protein. MKALLLLVLPWLSPANYIDNVGNLHFLYSELCKGASHYGLTKDRKRRSQDGCP DGCASLTATAPSPEVSAAATISLMTDEPGLDNPAYVSSAEDGQPAISPVDSGR SNRTRARPFERSTIRSRSFKKINRALSVLRRTKSGSAVANHADQGRESENTT APEVFPRLYHLIPDGEITSIKINRVDPSESLSIRLVGGSETPLVHIIIQHIYRDGVIA RDGRLLPGDIILKVNGMDISNVPHNYAVRLLRQPCQVLWLTVMREQKFRSRNN GQAPDAYRPRDDSFHVILNKSSPEEQLGIKLVKRVDEPGVFIFNVLDGGVAYRH GQLEENDRVLAINGHDLRYGSPESAHLIQASERRVHLVVSQRQVRQRSPDIFQ EAGWNSNGSWSPGPGERSNTPKPLHPTITCHEKVNIQKDPGESLGMTVAGG ASHREWDLPYVISVEPGGVISRDGRIKTGDILLNVDGVELTEVSRSEAVALLKR TSSSIVLKALEVKEYEPQEDCSSPAALDSNHNMAPPSDWSPSWVMWLELPRC LYNCKDIVLRRNTAGSLGFCIVGGYEEYNGNKPFFIKSIVEGTPAYNDGRIRCGD ILLAVNGRSTSGMIHACLARLLKELKGRITLTIVSWPGTFL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, Functional
Recommended Dilutions	Western Blot, Functional
Application Notes	This antibody is useful for Western Blot, Functional



Images

Western Blot: Lnx1 Antibody [H00084708-B01P] - Analysis of LNX1 expression in transfected 293T cell line by LNX1 polyclonal antibody. Lane 1: LNX1 transfected lysate(69.52 KDa). Lane 2: Non-transfected lysate.



Publications

Nayak D, Sivaraman J., et al. Structure of LNX1:Ubc13~Ubiquitin Complex Reveals the Role of Additional Motifs for the E3 Ligase Activity of LNX1. J Mol Biol 2018-02-26 [PMID: 29496391]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00084708-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

