

# Product Datasheet

## OSR1 Antibody (1D5) H00130497-M22

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00130497-M22](http://www.novusbio.com/H00130497-M22)

Updated 4/25/2019 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00130497-M22](http://www.novusbio.com/reviews/destination/H00130497-M22)



**H00130497-M22**

OSR1 Antibody (1D5)

| Product Information         |                                                                                                                                    |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 mg                                                                                                                             |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                   |
| Storage                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                       |
| Clonality                   | Monoclonal                                                                                                                         |
| Clone                       | 1D5                                                                                                                                |
| Preservative                | No Preservative                                                                                                                    |
| Isotype                     | IgG2b Kappa                                                                                                                        |
| Purity                      | IgG purified                                                                                                                       |
| Buffer                      | PBS (pH 7.4)                                                                                                                       |
| Product Description         |                                                                                                                                    |
| Host                        | Mouse                                                                                                                              |
| Gene ID                     | 130497                                                                                                                             |
| Gene Symbol                 | OSR1                                                                                                                               |
| Species                     | Human                                                                                                                              |
| Reactivity Notes            | Human. Other species not tested.                                                                                                   |
| Specificity/Sensitivity     | OSR1 (1D5)                                                                                                                         |
| Immunogen                   | OSR1 (NP_660303, 1 a.a. ~ 101 a.a) partial recombinant protein with GST tag.                                                       |
| Notes                       | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                 |
| Product Application Details |                                                                                                                                    |
| Applications                | Western Blot, ELISA                                                                                                                |
| Recommended Dilutions       | Western Blot 1:500, ELISA                                                                                                          |
| Application Notes           | Antibody Reactive against Recombinant Protein with GST tag on ELISA and Western Blot. GST tag alone is used as a negative control. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00130497-M22](http://www.novusbio.com/reviews/submit/H00130497-M22)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

