

# Product Datasheet

## PIWIL4 Antibody H00143689-B02P-50ug

Unit Size: 50 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00143689-B02P-50ug](http://www.novusbio.com/H00143689-B02P-50ug)

Updated 1/12/2021 v.20.1

### Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00143689-B02P-50ug](http://www.novusbio.com/reviews/destination/H00143689-B02P-50ug)




**H00143689-B02P-50ug**

## PIWIL4 Antibody

Product Information	
Unit Size	50 ug
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	143689
Gene Symbol	PIWIL4
Species	Human
Immunogen	PIWIL4 (ABM85351.1, 1 a.a. - 852 a.a.) full-length human protein. MSGRARVKARGIARSPSATEVGRIQASPLPRSVDLNSNEASSSNGFLGTSRIST NDKYGISSGDAGSTFMERGVKNKQDFMDLSICTREKLAHVRNCKTGSSGIPVK LVTNLFNLDFPQDWQLYQYHVITYIPDLASRRRLRIALLYSHSELSNKAKAFDGA ILFLSQKLEEKVTELSSETQRGETIKMTITLKRELPSSSPVCIQVFNIIFRKILKLSM YQIGRNFYNPSEPMEIPQHKLSLWPGFAISVSYFERKLLFSADVSYKVLNRNETV LEFMTALCQRTGLSCFTQTCEKQLIGLIVLTRYNNRTYSIDDIDWSVKPHTHFQK RDGTEITYVDYYKQQYDITVSDLNQPMLVSLKRRNDNSEAQLAHLIPELCFLT GLTDQATSDFQLMKAVAEKTRLSPSGRQRLARLVDNIQRNTNARFELETWGL HFGSQISLTGRIVPSEKILMQDHICQPVSAAADWSKDIRTCKILNAQSLNTWLILCS DRTEYVAESFLNCLRRVAGSMGFNVDPYPKIIVQENPAAFVRAIQQYVDPDVQL VMCILPSNQKTYYSIKKYLSSDCPVPSQCVLARTLNKQGMMSIATKIAMQMT CKLGGELWAVEIPLKSLMVVGDVCKDALSKDVMVVGCVASVNPRITRWFSRCI LQRTMTDVADCLKVFMTGALNKWYKYNHDLPARIIVYRAGVGDGQLKTLIEYEV PQLLSSVAESSSNTSSRLSVIVVRKKCMPRFFTEMNRTVQNPPLGTVVVDSEAT RNEWYDFYLISQVACRGTVSPTYYNVIYDDNGLKPDHMQRLLTFKLCHLYYNWP GIVSVPAPCQYAHKLTFLVAQSIHKEPSLELANHLFYL
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot

## Images

Western Blot: PIWIL4 Antibody [H00143689-B02P-50ug] - Western Blot analysis of PIWIL4 expression in transfected 293T cell line by PIWIL4 MaxPab polyclonal antibody. Lane 1: PIWIL4 transfected lysate(93.72 kDa). Lane 2: Non-transfected lysate.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00143689-B02P-50ug](http://www.novusbio.com/reviews/submit/H00143689-B02P-50ug)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

