

Product Datasheet

ZNF342 Antibody H00162979-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00162979-B01P

Updated 1/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00162979-B01P



H00162979-B01P**ZNF342 Antibody**

Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein G purified
Buffer	PBS (pH 7.4)


Product Description	
Host	Mouse
Gene ID	162979
Gene Symbol	ZNF296
Species	Human
Specificity/Sensitivity	This product is specific for Human ZNF296 purified MaxPab mouse polyclonal antibody (B01P) [Gene ID: 162979].
Immunogen	ZNF296 (NP_660331.1, 1 a.a. - 475 a.a.) full-length human protein. MSRRKAGSAPRRVEPAPAANPDDEMEMQDLVIELKPEPDAQPQQAPRLGPFS PKEVSSAGRFGGEPHHSPGPMPAGAALLALGPRNPWTLWTPLTPNYPDRQP WTDKHPDLLTCGRCLQTFPLEAITAFMDHKKLGQCQLFRGPSRGQGSEREELKA LSCLRCGKQFTVAWKLLRHAQWDHGLSIYQTESEAPEAPLLGLAEVAAAVSAV VGPAAEAKSPRASGSLTRRSPTCPVCKKTLSSFSNLKVHMRSHTERPYAC DQCPYACAQSSKLNHRHKKTHRQVPPQSPLMADTSQEQASAAPPEPAVHAAAP TSTLPCSGGEGAGAAATAGVQEPGAPGSGAQAGPGGDTWGAIITEQRTDPA NSQKASPKKMPKSGGKSRGPGGSCEFCGKHFTNSSNLTVHRRSHTGERPYT CEFCNYACAQSSKLNHRHRMHGMTPGSTRFECPHCHVPFGLRATLDKHLRQK HPEAAGEA
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin
Application Notes	Mouse polyclonal antibody raised against a full-length human ZNF296 protein.




Images


Western Blot: ZNF342 Antibody [H00162979-B01P] - Analysis of ZNF296 expression in Hela S3 NE.




Immunocytochemistry/Immunofluorescence: ZNF342 Antibody [H00162979-B01P] - Analysis of purified antibody to ZNF342 on HeLa cell. (antibody concentration 10 ug/ml)



Immunohistochemistry-Paraffin: ZNF342 Antibody [H00162979-B01P] - Analysis of purified antibody to ZNF342 on formalin-fixed paraffin-embedded human lymph node. (antibody concentration 3 ug/ml)



Western Blot: ZNF342 Antibody [H00162979-B01P] - Analysis of ZNF296 expression in transfected 293T cell line by ZNF296 polyclonal antibody. Lane 1: ZNF342 transfected lysate(52.25 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis.
Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00162979-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

