

# Product Datasheet

## HIPK1 Antibody (1D6) - Azide and BSA Free H00204851-M06

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00204851-M06](http://www.novusbio.com/H00204851-M06)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00204851-M06](http://www.novusbio.com/reviews/destination/H00204851-M06)



**H00204851-M06**

HIPK1 Antibody (1D6) - Azide and BSA Free

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 mg                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| Clonality           | Monoclonal                                                                                                       |
| Clone               | 1D6                                                                                                              |
| Preservative        | No Preservative                                                                                                  |
| Isotype             | IgG2a Kappa                                                                                                      |
| Purity              | IgG purified                                                                                                     |
| Buffer              | In 1x PBS, pH 7.4                                                                                                |


| Product Description     |                                                                                                                                                                                                                                 |
|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description             | Novus Biologicals Mouse HIPK1 Antibody (1D6) - Azide and BSA Free (H00204851-M06) is a monoclonal antibody validated for use in IHC and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.              |
| Host                    | Mouse                                                                                                                                                                                                                           |
| Gene ID                 | 204851                                                                                                                                                                                                                          |
| Gene Symbol             | HIPK1                                                                                                                                                                                                                           |
| Species                 | Human                                                                                                                                                                                                                           |
| Reactivity Notes        | Human. Other species not tested.                                                                                                                                                                                                |
| Specificity/Sensitivity | HIPK1 (1D6)                                                                                                                                                                                                                     |
| Immunogen               | HIPK1 (AAH28408, 330 a.a. ~ 430 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.<br>CTPLMVATLHPQVATITPQYAVPFTLSAAGRPAALVEQTA AVLQAWPGGTQQIL<br>LPSTWQQQLPGVALHNSVQPTAMIPAMGSGGQLADWRNAHSHGNQYS |
| Notes                   | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                              |

| Product Application Details |                                                                                                                 |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------|
| Applications                | Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry, Sandwich ELISA                                      |
| Recommended Dilutions       | ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin, Sandwich ELISA                                      |
| Application Notes           | Antibody reactive against recombinant protein on ELISA. Has also been used for immunohistochemistry (paraffin). |




## Images

Immunohistochemistry-Paraffin: HIPK1 Antibody (1D6) [H00204851-M06] - Analysis of monoclonal antibody to HIPK1 on formalin-fixed paraffin-embedded human salivary gland. Antibody concentration 3 ug/ml



Sandwich ELISA: HIPK1 Antibody (1D6) [H00204851-M06] - Detection limit for recombinant GST tagged HIPK1 is 3 ng/ml as a capture antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00204851-M06**

---

|                  |                                                         |
|------------------|---------------------------------------------------------|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-96981-0.5mg | Mouse IgG2a Kappa Isotype Control (M2AK)                |
| NBP1-82981PEP    | HIPK1 Recombinant Protein Antigen                       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00204851-M06](http://www.novusbio.com/reviews/submit/H00204851-M06)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

