

Product Datasheet

HEATR9 Antibody H00256957-B01P

Unit Size: 0.05 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00256957-B01P

Updated 7/19/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00256957-B01P




H00256957-B01P

HEATR9 Antibody


Product Information	
Unit Size	0.05 mg
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	256957
Gene Symbol	C17ORF66
Species	Human
Specificity/Sensitivity	Reacts with chromosome 17 open reading frame 66.
Immunogen	C17orf66 (AAH33734, 1 a.a. - 219 a.a.) full-length human protein. MAYEKSTDISDVSRSMFLYPWLEYPDKTKELRKAMAPVHPLPLSCYQMPKEEFP PSPECWRQHPSKPNVSPYCYFKKPEIYTHWHDLYDQREEREA EKMLRKMRD DCRYIKEVHQTHIKMFHLPMSKLTIKSEMRSRPLEPTQDPLKWQRLRELTKSLE SPREDEQFYAAQALGCLRISDKFVMEALQQVVSLRTAGVQERRKQGLRCVAC GHRRGSPGC
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/Immunofluorescence
Application Notes	This antibody is reactive against tissue and transfected lysate in western blot, immunofluorescence, and as a detection antibody in ELISA.

Images


Western Blot: C17orf66 Antibody [H00256957-B01P] - Analysis of C17orf66 expression in human kidney.



Immunocytochemistry/Immunofluorescence: C17orf66 Antibody [H00256957-B01P] - Analysis of purified antibody to C17orf66 on HeLa cell. (antibody concentration 10 ug/ml)



Western Blot: C17orf66 Antibody [H00256957-B01P] - Analysis of C17orf66 expression in transfected 293T cell line by C17orf66 polyclonal antibody. Lane 1: C17orf66 transfected lysate(24.09 KDa). Lane 2: Non-transfected lysate.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00256957-B01P

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97019-5mg	Mouse IgG Isotype Control
H00256957-P01-10ug	Recombinant Human HEATR9 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00256957-B01P

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

