

# Product Datasheet

## TTC9C Antibody (7H10) - Azide and BSA Free H00283237-M06-100ug

Unit Size: 100 ug

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00283237-M06-100ug](http://www.novusbio.com/H00283237-M06-100ug)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00283237-M06-100ug](http://www.novusbio.com/reviews/destination/H00283237-M06-100ug)



**H00283237-M06-100ug**

TTC9C Antibody (7H10) - Azide and BSA Free

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 100 ug                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| Clonality           | Monoclonal                                                                                                       |
| Clone               | 7H10                                                                                                             |
| Preservative        | No Preservative                                                                                                  |
| Isotype             | IgG2a Kappa                                                                                                      |
| Purity              | Protein A or G purified                                                                                          |
| Buffer              | In 1x PBS, pH 7.4                                                                                                |


| Product Description |                                                                                                                                                                                                                                                                                                                     |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Mouse TTC9C Antibody (7H10) - Azide and BSA Free (H00283237-M06-100ug) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                            |
| Host                | Mouse                                                                                                                                                                                                                                                                                                               |
| Gene ID             | 283237                                                                                                                                                                                                                                                                                                              |
| Gene Symbol         | TTC9C                                                                                                                                                                                                                                                                                                               |
| Species             | Human                                                                                                                                                                                                                                                                                                               |
| Immunogen           | TTC9C (NP_776171.1, 1 a.a. ~ 171 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.<br>MEKRLQEAQLYKEEGNQRYREGKYRDAVSRYHRALLQLRGLDPSLPSPPLPNL<br>GPQGPALTPEQENILHTTQTDCCYNNLAACLLQMEPVNYERVREYSQKVLERQP<br>DNAKALYRAGVAFFHLQDYDQARHYLLAAVNRQPKDANVRRYLQLTQSELSS<br>YHRKEKQLYLGMFG |
| Notes               | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                  |

| Product Application Details |                     |
|-----------------------------|---------------------|
| Applications                | Western Blot, ELISA |
| Recommended Dilutions       | Western Blot, ELISA |




## Images

Western Blot: TTC9C Antibody (7H10) [H00283237-M06-100ug] - Detection against Immunogen (46.4 KDa) .



ELISA: TTC9C Antibody (7H10) [H00283237-M06-100ug] - Detection limit for recombinant GST tagged TTC9C is 1 ng/ml as a capture antibody.



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to H00283237-M06-100ug**

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| HAF007             | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B            | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-96981-0.5mg   | Mouse IgG2a Kappa Isotype Control (M2AK)                |
| H00283237-P01-10ug | Recombinant Human TTC9C GST (N-Term) Protein            |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00283237-M06-100ug](http://www.novusbio.com/reviews/submit/H00283237-M06-100ug)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

