

Product Datasheet

Recombinant Human TMPRSS7 Protein H00344805-G01-10ug

Unit Size: 10 ug

Store at -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/H00344805-G01-10ug

Updated 11/4/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/H00344805-G01-10ug



H00344805-G01-10ug

Recombinant Human TMPRSS7 Protein

Product Information	
Unit Size	10 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -80C. Avoid freeze-thaw cycles.
Preservative	Glycerol
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	25 mM Tris-HCl pH8.0 in 2% glycerol.
Target Molecular Weight	63 kDa
Product Description	
Description	<p>An un-tageed recombinant protein corresponding to full length human TMPRSS7 ORF (ADR83226.1)</p> <p>Source: <i>Wheat Germ (in vitro) with proprietary liposome technology</i></p> <p>Amino Acid Sequence: MCHFKLVAIVGYLIRLSIKSIEADNCVTDSLTIYDSELLPIRSSILYRICEPTRTLM SFVSTNNLMLLVTFKSPHIRRLSGIRAYFEVIPEQKCENTVLVKDITGFEGKISSPY YPSYPPKCKCTWKFQTSLSLTLGIALKFYNYSITKKSMKGCEHGWWEINEHMY CGSYMDHQITFRVPSPLVHIQLQCSSLSDKPLLAEYGSYNISQPCPVGSFRCS SGLCVPQAQRCDGVNDCFDESDELFCVSPQPACNTSSFRQHGPLICDGFDC ENGRDEQNCTQSIPCNNRTFKCGNDICFRKQNAKCDGTVDPCDGSDEEGCTC SRSSALHRIIGGTDLEGGWPWQVSLHFVGSAYCGASVISREWLLSAAHCFH GNRLSDPTPWAHLGMYVQGNKAFVSPVRRIVVHEYYNSQTFDYDIALQLSIA WPETLKQLIQPICIPPTGQRVRSGEKCVWVGWRRHEADNKGSLVLQQAEEV LIDQTLCVSTYGIITSRMLCAGIMSGKRDACKGDSGGPLSCRKSDGKWILTGV SWGHCGRPNFPGVYTRVSNFVPIHKYVPSLL</p>
Gene ID	344805
Gene Symbol	TMPRSS7
Species	Human
Preparation Method	in vitro wheat germ expression system with proprietary liposome technology
Details of Functionality	This protein was produced in an in vitro wheat germ expression system that should preserve correct conformational folding that is necessary for biological function. While it is possible that this protein could display some level of activity, the functionality of this protein has not been explicitly measured or validated.
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Product Application Details	
Applications	Immunoaffinity Purification
Recommended Dilutions	Immunoaffinity Purification



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to H00344805-G01-10ug

NBP1-90893PEP	TMPRSS7 Recombinant Protein Antigen
NB600-41532	Albumin Antibody
NBP1-90893	TMPRSS7 Antibody
AF3946	Matriptase/ST14 Antibody

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/H00344805-G01-10ug

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

