

# Product Datasheet

## Surfactant Protein A Antibody Pair H00653509-AP43

Unit Size: 1 Set

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00653509-AP43](http://www.novusbio.com/H00653509-AP43)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00653509-AP43](http://www.novusbio.com/reviews/destination/H00653509-AP43)



**H00653509-AP43****Surfactant Protein A Antibody Pair****Product Information**

<b>Unit Size</b>	1 Set
<b>Concentration</b>	Concentration of individual antibodies may be found on the vial label. If unlisted please contact technical services.
<b>Storage</b>	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
<b>Preservative</b>	No Preservative
<b>Purity</b>	Immunogen affinity purified
<b>Buffer</b>	1x PBS, pH 7.4

**Product Description**


<b>Description</b>	Quality control test: Standard curve using recombinant protein (H00653509-P01) as an analyte. This antibody pair set comes with matched antibody pair to detect and quantify protein level of human SFTPA1.
<b>Gene ID</b>	653509
<b>Gene Symbol</b>	SFTPA1
<b>Species</b>	Human
<b>Kit Components</b>	Capture antibody: mouse monoclonal anti-SFTPA1 (100 ug), Detection antibody: biotinylated mouse monoclonal anti-SFTPA1 IgG1 Kappa (50 ug)
<b>Notes</b>	This product is produced by and distributed for Abnova, a company based in Taiwan.
<b>Standard Curve Range</b>	3 ng/ml to 100 ng/ml
<b>Assay Type</b>	Sandwich ELISA

**Product Application Details**

<b>Applications</b>	Sandwich ELISA
<b>Recommended Dilutions</b>	Sandwich ELISA 3 ng/ml - 100 ng/ml

**Images**

Sandwich ELISA: Surfactant Protein A Antibody Pair [H00653509-AP43]  
- Detection sensitivity ranging from 3 ng/ml to 100 ng/ml.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to H00653509-AP43**

---

H00653509-P01-10ug	Recombinant Human Surfactant Protein A GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
NBP2-12928	Surfactant Protein A Antibody
M6000B-1	IL-6 [HRP]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Antibody Pairs are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00653509-AP43](http://www.novusbio.com/reviews/submit/H00653509-AP43)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

