

Product Datasheet

Human Integrin alpha 2 beta/Integrin beta 3 ELISA Kit (Colorimetric) KA0477

Unit Size: 1 Kit

Storage of components varies. See protocol for specific instructions.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/KA0477

Updated 4/1/2018 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/KA0477



KA0477**Human Integrin alpha 2 beta/Integrin beta 3 ELISA Kit (Colorimetric)**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Storage of components varies. See protocol for specific instructions.
Product Description	
Species	Human
Reactivity Notes	Human
Kit Components	Human GPIIb/IIIa Microplate: A 96-well polystyrene microplate coated with a polyclonal antibody against GPIIb/IIIa., Sealing Tapes: Pre-cut, pressure-sensitive sealing tapes that can be cut to fit the format of the individual assay., Human GPIIb/IIIa Standard: Human platelet GPIIb/IIIa in a buffered protein base (lyophilized)., Biotinylated Human GPIIb/IIIa Antibody (100x): A 100-fold concentrated biotinylated polyclonal antibody against human GPIIb/IIIa., MIX Diluent Concentrate (10x): A 10-fold concentrated buffered protein base., Wash Buffer Concentrate (20x): A 20-fold concentrated buffered surfactant., Streptavidin-Peroxidase Conjugate (SP Conjugate): A 100-fold concentrate., Chromogen Substrate: A ready-to-use stabilized peroxidase chromogen substrate tetramethylbenzidine., Stop Solution: A 0.5 N hydrochloric acid to stop the chromogen substrate reaction.
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Standard Curve Range	1.25 to 40 ng/mL
Sensitivity	1 ng/mL
Inter Assay	CV% < 7.70%
Intra Assay	CV% < 4.70%
Assay Type	Sandwich ELISA
Spiking Recovery	99.00%
Suitable Sample Type	Cell Lysate, Plasma, Platelet
Sample Volume	50 uL
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA
Application Notes	ITGA2B/ITGB3 (Human) ELISA Kit is a sandwich enzyme immunoassay for the quantitative measurement of human ITGA2B/ITGB3.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/KA0477

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

