

Product Datasheet

ANCA screen ELISA Kit (Colorimetric) KA1083

Unit Size: 1 Kit

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/KA1083

Updated 6/27/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/KA1083



KA1083**ANCA screen ELISA Kit (Colorimetric)**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C.
Product Description	
Description	Quality control test: Standard curve. ANCAbScreen ELISA Kit is a assay intended to screen for the presence of IgG class autoantibodies against PR3 and MPO in human serum or plasma.
Kit Components	Qty.1: Divisible microplate consisting of 12 modules of 8 wells each, coated with a mixture of highly purified antigens, PR3 and MPO. Ready to use., 3 vials, 1.5 ml each: Anti-ANCA controls in a serum/buffer matrix (PBS, NaN ₃ <0.1% (w/w)). Negative Control (A), Cut-Off Control (B) and Positive Control (C). Ready to use., 1 vial, 20 ml: Sample buffer (Tris, NaN ₃ <0.1% (w/w)), yellow, concentrate (5x), 1 vial, 15 ml: Enzyme conjugate solution (PBS, PROCLIN 300 <0.5% (v/v)), (light red) containing polyclonal rabbit anti-human IgG; labelled with horseradish peroxidase. Ready to use., 1 vial, 15 ml: TMB substrate solution. Ready to use., 1 vial, 15 ml: Stop solution (1 M hydrochloric acid). Ready to use., 1 vial, 20 ml: Wash solution (PBS, NaN ₃ <0.1% (w/w)), concentrate (50x)
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Sensitivity	not applicable
Inter Assay	CV% < 6.3-9.8%
Intra Assay	CV% < 6.3-9.1%
Suitable Sample Type	Plasma, Serum
Sample Volume	100 uL (prediluted samples)
Product Application Details	
Applications	ELISA
Recommended Dilutions	ELISA
Application Notes	ANCAbScreen ELISA Kit is a qualitative assay intended to screen for the presence of IgG class autoantibodies against PR3 and MPO in human serum or plasma.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/KA1083

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

