

Product Datasheet

Human Rheumatoid Factor ELISA Kit (Colorimetric) KA1289

Unit Size: 1 Kit

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/KA1289

Updated 4/12/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/KA1289



KA1289**Human Rheumatoid Factor ELISA Kit (Colorimetric)**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C.
Product Description	
Description	Rheumatoid factor IgA ELISA Kit is indirect solid phase enzyme immunoassay for the quantitative measurement of rheumatoid factor IgA in serum and plasma.
Species	Human
Kit Components	Divisible microplate consisting of 12 modules of 8 wells each. Ready to use., Calibrators A-E (0, 15, 50, 150, 500 U/mL), containing rheumatoid factor in a serum/buffer matrix (PBS, BSA, detergent, NaN ₃ 0.09%), yellow. Ready to use., Control positive (1) and negative (2), containing rheumatoid factor in a serum/buffer matrix (PBS, BSA, detergent, NaN ₃ 0.09%), yellow. Ready to use. The concentration is specified on the certificate of analysis., Sample buffer P, containing PBS, BSA, detergent, preservative NaN ₃ 0.09%, yellow, concentrate (5x), Enzyme conjugate, light red, containing anti-human IgA antibodies, HRP labelled; PBS, BSA, detergent, preservative Proclin 300 0.05%, Ready to use., TMB substrate, containing 3,3',5,5'- Tetramethylbenzidine, colorless. Ready to use., Stop solution, contains acid. Ready to use., Wash solution, containing Tris, detergent, preservative NaN ₃ 0.09%, 50x conc.
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Standard Curve Range	0-500 U/mL
Sensitivity	1 U/mL
Inter Assay	CV% < 9.10%
Intra Assay	CV% < 6.20%
Assay Type	Indirect ELISA
Spiking Recovery	80-113%
Suitable Sample Type	Plasma, Serum
Sample Volume	100 uL (prediluted)
Product Application Details	
Applications	ELISA, Quantification
Recommended Dilutions	ELISA, Quantification
Application Notes	Useful for quantification.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/KA1289

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

