

Product Datasheet

Metanephrine ELISA Kit (Colorimetric) KA1889

Unit Size: 1 Kit

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/KA1889

Updated 12/30/2020 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/KA1889



KA1889**Metanephrine ELISA Kit (Colorimetric)**


| Product Information | |
|----------------------------|----------------------------------------------------------------------------------------------------------|
| Unit Size | 1 Kit |
| Concentration | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| Storage | Store at 4C. |

| Product Description | |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description | Metanephrine ELISA Kit is an enzyme immunoassay for the quantitative determination of free Metanephrine in plasma. |
| Kit Components | Adhesive Foil, Wash Buffer Concentrate (50X), Enzyme Conjugate, Substrate, Stop Solution, Metanephrine Microtiter Strips, Metanephrine Antiserum, Adjustment Buffer, Assay Buffer, Acylation Concentrate, Extraction Plate, Cleaning Concentrate (25x), Elution Buffer, Equalizing-Reagent |
| Notes | This product is produced by and distributed for Abnova, a company based in Taiwan. |
| Standard Curve Range | 15.1 to 3600 pg/mL |
| Sensitivity | 14.9 pg/mL |
| Inter Assay | CV% < 9.47% |
| Intra Assay | CV% < 9.77% |
| Assay Type | Quantification |
| Spiking Recovery | 1.129 |
| Suitable Sample Type | Plasma |
| Sample Volume | 50 uL (extracted samples) |

| Product Application Details | |
|------------------------------------|----------------------------|
| Applications | ELISA, Quantification |
| Recommended Dilutions | ELISA, Quantification |
| Application Notes | Useful for quantification. |

Images

ELISA: Metanephrine ELISA Kit (Colorimetric) [KA1889] - These standard curves are provided for demonstration only. A standard curve should be generated for each set of samples assayed.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/KA1889

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

