

# Product Datasheet

## Porcine Plasminogen ELISA Kit (Colorimetric) KA1959

Unit Size: 1 Kit

Store at 4C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/KA1959](http://www.novusbio.com/KA1959)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/KA1959](http://www.novusbio.com/reviews/destination/KA1959)



**KA1959****Porcine Plasminogen ELISA Kit (Colorimetric)**

| <b>Product Information</b>         |  |
|------------------------------------|--|
| <b>Unit Size</b>                   | 1 Kit  |
| <b>Concentration</b>               | Concentration is not relevant for this product. Please see the protocols for proper use of this product.   |
| <b>Storage</b>                     | Store at 4C.   |
| <b>Product Description</b>         |  |
| <b>Description</b>                 | PLG (Pig) ELISA Kit is a highly sensitive two-site enzyme-linked immunoassay for the quantitative determination of PLG in pig biological fluids. |
| <b>Gene ID</b>                     | 5340   |
| <b>Gene Symbol</b>                 | PLG  |
| <b>Species</b>                     | Porcine  |
| <b>Notes</b>                       | This product is produced by and distributed for Abnova, a company based in Taiwan.   |
| <b>Standard Curve Range</b>        | 6.25 to 200 ng/mL  |
| <b>Sensitivity</b>                 | 1.142  |
| <b>Assay Type</b>                  | Two-site enzyme-linked immunoassay   |
| <b>Suitable Sample Type</b>        | Biological Sample  |
| <b>Sample Volume</b>               | 100 uL   |
| <b>Product Application Details</b> |  |
| <b>Applications</b>                | ELISA, Quantification  |
| <b>Recommended Dilutions</b>       | ELISA, Quantification  |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to KA1959**

---

|                 |   |
|-----------------|---|
| NBP1-86015PEP   | Plasminogen Recombinant Protein Antigen |
| 210-TA-005      | TNF-alpha [Unconjugated]                |
| NB600-930-100ug | Plasminogen Antibody                    |
| M6000B-1        | IL-6 [HRP]                              |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/KA1959](http://www.novusbio.com/reviews/submit/KA1959)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

