

Product Datasheet

Mouse Transthyretin/Prealbumin ELISA Kit (Colorimetric) KA2070

Unit Size: 1 Kit

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/KA2070

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/KA2070



KA2070**Mouse Transthyretin/Prealbumin ELISA Kit (Colorimetric)**

Product Information	
Unit Size	1 Kit
Concentration	Concentration is not relevant for this product. Please see the protocols for proper use of this product.
Storage	Store at 4C.
Product Description	
Description	Prealbumin (Mouse) ELISA Kit is an enzyme immunoassay for the quantitative determination of Prealbumin in mouse biological samples.
Gene ID	7276
Gene Symbol	TTR
Species	Mouse
Kit Components	Diluent Concentrate (Running Buffer): One bottle containing 5X concentrated diluent running buffer., Wash Solution Concentrate: One bottle containing 20X concentrated wash solution., Enzyme-Antibody Conjugate 100X: One vial containing affinity-purified anti- PRE antibody conjugated with horseradish peroxidase in a stabilizing buffer., Chromogen-Substrate Solution: One vial containing 3,3',5,5'-tetramethylbenzidine (TMB) and hydrogen peroxide in citric acid buffer at pH 3.3., Stop Solution: One vial containing 0.3 M sulfuric acid. WARNING: Avoid contact with skin., ANTI-PRE ELISA MICRO PLATE: Twelve removable eight (8) well micro well strips in well holder frame. Each well is coated with affinity purified anti-PRE., MOUSE PRE CALIBRATOR: One vial containing a lyophilized Mouse PRE calibrator.
Notes	This product is produced by and distributed for Abnova, a company based in Taiwan.
Standard Curve Range	15.63 to 500 ng/mL
Sensitivity	3.224
Assay Type	Quantification
Suitable Sample Type	Biological samples
Sample Volume	100 uL
Product Application Details	
Applications	ELISA, Quantification
Recommended Dilutions	ELISA, Quantification



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to KA2070

NBP3-14788-100ug	Recombinant Human Transthyretin/Prealbumin Y78F Variant, Filaments Protein
7954-GM-010/CF	GM-CSF [Unconjugated]
NBP2-52575	Transthyretin/Prealbumin Antibody (2E10C5) - BSA Free
NBP2-13075	beta Amyloid Antibody (MOAB-2) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/KA2070

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

