

# Product Datasheet

## Protein Disulfide Isomerase/P4HB Antibody - Azide and BSA Free NB100-1921

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

**Reviews: 2 Publications: 6**

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-1921](http://www.novusbio.com/NB100-1921)

Updated 9/9/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-1921](http://www.novusbio.com/reviews/destination/NB100-1921)



**NB100-1921**


Protein Disulfide Isomerase/P4HB Antibody - Azide and BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Unpurified
Buffer	Whole antisera
Product Description	
Description	Novus Biologicals Rabbit Protein Disulfide Isomerase/P4HB Antibody - Azide and BSA Free (NB100-1921) is a polyclonal antibody validated for use in IHC, WB, ICC/IF and IP. Anti-Protein Disulfide Isomerase/P4HB Antibody: Cited in 6 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5034
Gene Symbol	P4HB
Species	Human, Mouse, Rat, Bovine, Primate
Marker	Endoplasmic Reticulum Marker
Immunogen	Purified full length P4hb isolated from rat liver.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Western Blot 1:10000-1:20000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:5000-1:20000, Immunoprecipitation 1-3 ug, Immunohistochemistry-Paraffin 1:100-1:500, Immunohistochemistry-Frozen 1:100-1:1000
Application Notes	This P4hb antibody is useful for Immunocytochemistry/Immunofluorescence, Immunohistochemistry frozen and paraffin embedded sections, Immunoprecipitation and Western Blot.




## Images


Normal human brain.



ICC/IF Confocal analysis of 3T3 cells using P4HB antibody (NB100-1921, 1:5). An Alexa Fluor 488-conjugated Goat to rabbit IgG was used as secondary antibody (green). Actin filaments were labeled with Alexa Fluor 568 phalloidin (red). DAPI was used to stain the cell nuclei (blue).



IHC analysis of a formalin fixed paraffin-embedded (FFPE) human brain using 1:250 conc. of P4HB antibody on a Bond Rx autostainer (Leica Biosystems). The assay involved 20 minutes of heat induced antigen retrieval (HIER) using 10mM sodium citrate buffer (pH 6.0) and endogenous peroxidase quenching with peroxide block. The sections were incubated with primary antibody for 30 minutes and Bond Polymer Refine Detection (Leica Biosystems) with DAB was used for signal development followed by counterstaining with hematoxylin. Whole slide scanning and capturing of representative images was performed using Aperio AT2 (Leica Biosystems). Cytoplasmic staining of P4HB in neurons was observed. Staining was performed by Histowiz.



P4HB staining in Hela cells detected with a Dylight 488 labeled secondary antibody (Green) with Hoechst 33258 nuclear counterstain (Blue).



## Publications

Oleinik N, Kim J, Roth BM et al. Mitochondrial protein import is regulated by p17/PERMIT to mediate lipid metabolism and cellular stress Sci Adv 2019-09-01 [PMID: 31535025] (WB, Human)

Gamberucci A, Marcolongo P, Nemeth CE et al. GLUT10-Lacking in Arterial Tortuosity Syndrome-Is Localized to the Endoplasmic Reticulum of Human Fibroblasts Int J Mol Sci 2017-08-22 [PMID: 28829359] (ICC/IF, Human)

Omachi K, Kamura M, Teramoto K et al. A Split-Luciferase-Based Trimer Formation Assay as a High-throughput Screening Platform for Therapeutics in Alport Syndrome. Cell Chem Biol. 2018-02-21 [PMID: 29526710] (WB)

Jun HS, Lee YM, Song K et al. G-CSF improves murine G6PC3-deficient neutrophil function by modulating apoptosis and energy homeostasis Blood 2011-04-07 [PMID: 21292774]

Azad AK, Torrelles JB, Schlesinger LS. Mutation in the DC-SIGN cytoplasmic triacidic cluster motif markedly attenuates receptor activity for phagocytosis and endocytosis of mannose-containing ligands by human myeloid cells. J Leukoc Biol;84(6):1594-603. 2008-12-01 [PMID: 18772280] (ICC/IF, Human)

Safran, M et al. Thyroid hormone-dependent redistribution of the 55-kilodalton monomer of protein disulfide isomerase in cultured glial cells. Endocrinology;131(5):2413-8. 1992-11-01 [PMID: 1425439]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-1921**

---

NB820-60576	Human Brain Frontal Lobe Membrane Tissue Lysate (Adult Membrane Normal)
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-1921](http://www.novusbio.com/reviews/submit/NB100-1921)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

