

Product Datasheet

SOD2/Mn-SOD Antibody - BSA Free NB100-1969

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 5

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-1969

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-1969



NB100-1969

SOD2/Mn-SOD Antibody - BSA Free

Product Information	
Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.4), 50% Glycerol


Product Description	
Description	Novus Biologicals Rabbit SOD2/Mn-SOD Antibody - BSA Free (NB100-1969) is a polyclonal antibody validated for use in IHC, WB, ELISA and IP. Anti-SOD2/Mn-SOD Antibody: Cited in 5 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	6648
Gene Symbol	SOD2
Species	Human, Mouse, Rat, Porcine, Bovine, Canine, Chicken, Guinea Pig, Hamster, Gerbil, Monkey, Rabbit, Sheep, Xenopus
Reactivity Notes	Mouse reactivity reported in scientific literature (PMID: 24578415). Immunogen displays sequence identity for non-tested species: Rat, Bovine, Canine, Chicken, Gerbil, Guinea pig, Porcine, Hamster, Monkey, Rabbit, Sheep, Xenopus, Frog and Wood Frog (<i>Rana sylvatica</i>).
Specificity/Sensitivity	Detects approx 25 kDa.
Immunogen	Recombinant Human Mn SOD Protein

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation
Recommended Dilutions	Western Blot 1:5000, ELISA 1:100 - 1:2000, Immunohistochemistry 1:100, Immunoprecipitation 1:10 - 1:500, Immunohistochemistry-Paraffin 1:100, Immunohistochemistry-Frozen 1:100
Application Notes	0.2 ug/mL of SOD (Mn) Antibody was sufficient for detection of Mn SOD in 20 ug of rat brain tissue extract by colorimetric immunoblot analysis using Goat anti-mouse IgG:AP as the secondary Antibody. Use in Immunohistochemistry-Frozen reported in scientific literature (PMID 24578415).



Images

Western Blot: SOD2/Mn-SOD Antibody [NB100-1969] - Western blot analysis of Mouse Brain, Rat Brain, and Rat Brain Membrane showing detection of 24.7kDa SOD (Mn) protein using Rabbit Anti-SOD (Mn) Polyclonal Antibody (NB100-1969). Lane 1: MW Ladder. Lane 2: Mouse Brain (15 ug). Lane 3: Rat Brain (15 ug). Lane 4: Rat Brain Membrane (15 ug). Block: 5% Skim Milk powder in TBST. Primary Antibody: Rabbit Anti-SOD (Mn) Polyclonal Antibody (NB100-1969) at 1:1000 for O/N at 4C. Secondary Antibody: Goat anti-rabbit IgG:HRP at 1:5000 for 1 hour at RT with shaking. Color Development: Chemiluminescent for HRP (Moss) for 5 min in RT. Predicted/Observed Size: 24.7kDa.



Publications

Lakhani R, Vogel KR, Till A et al. Defects in GABA metabolism affect selective autophagy pathways and are alleviated by mTOR inhibition. *EMBO Mol Med* 2014-02-28 [PMID: 24578415] (IHC-Fr, Mouse)

Kalash R, Epperly MW, Goff J et al. Amelioration of Radiation-Induced Pulmonary Fibrosis by a Water-Soluble Bifunctional Sulfoxide Radiation Mitigator (MMS350). *Radiat Res*. 2013-11-01 [PMID: 24125487]

Epperly MW, Chaillet JR, Kalash R et al. Conditional Radioresistance of tet-Inducible Manganese Superoxide Dismutase Bone Marrow Stromal Cell Lines. *Radiat Res* 2013-08-01 [PMID: 23862693]

Dumont M, Stack C, Elipenahli C et al. Behavioral deficit, oxidative stress, and mitochondrial dysfunction precede tau pathology in P301S transgenic mice. *FASEB J*. 2011-08-08 [PMID: 21825035]

Lee DW, Rajagopalan S, Siddiq A et al. Inhibition of Prolyl Hydroxylase Protects against 1-Methyl-4-phenyl-1,2,3,6-tetrahydropyridine-induced Neurotoxicity: MODEL FOR THE POTENTIAL INVOLVEMENT OF THE HYPOXIA-INDUCIBLE FACTOR PATHWAY IN PARKINSON DISEASE. *J Biol Chem*;284(42):29065-29076. 2009-01-01 [PMID: 19679656]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-1969

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-52079-0.1mg	Recombinant Mouse SOD2/Mn-SOD His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-1969

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

