

Product Datasheet

Dynamin Antibody (D5) - Azide and BSA Free NB100-1993

Unit Size: 0.05 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-1993

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-1993



NB100-1993


Dynamamin Antibody (D5) - Azide and BSA Free

Product Information	
Unit Size	0.05 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	D5
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS (pH 7.2), 0.1mM PMSF and 50% Glycerol
Product Description	
Description	Novus Biologicals Mouse Dynamamin Antibody (D5) - Azide and BSA Free (NB100-1993) is a monoclonal antibody validated for use in IHC, WB, Flow and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	1759
Gene Symbol	DNM1
Species	Human, Mouse, Rat, Bovine, Chicken
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information.
Immunogen	Recombinant His6-tagged human dynamamin.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 0.5 ug/mL, Flow Cytometry 10 ug/ml, Immunohistochemistry 1:50, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Paraffin 1:50
Application Notes	Detects a band of ~100 kDa by Western blot; recommended dilutions for colorimetric detection.




Images

Western Blot: Dynamin Antibody (D5) [NB100-1993] - Analysis of chicken SL29 cell Lysate



Western Blot: Dynamin Antibody (D5) [NB100-1993] - Western blot analysis of Dynamin, mAb (D5): Lane 1: MW marker, Lane 2: Mouse Brain, Lane 3: Rat Brain





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-1993

NB820-59657	Mouse Brain Whole Tissue Lysate (Adult Whole Normal)
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-1993

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

