

# Product Datasheet

## HIF-3 alpha Antibody [DyLight 650] NB100-2287C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-2287C](http://www.novusbio.com/NB100-2287C)

Updated 9/10/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-2287C](http://www.novusbio.com/reviews/destination/NB100-2287C)



**NB100-2287C**

HIF-3 alpha Antibody [DyLight 650]

| <b>Product Information</b>         |   |
|------------------------------------|---|
| <b>Unit Size</b>                   | 0.1 ml  |
| <b>Concentration</b>               | Please see the vial label for concentration. If unlisted please contact technical services.   |
| <b>Storage</b>                     | Store at 4C in the dark.  |
| <b>Clonality</b>                   | Polyclonal  |
| <b>Preservative</b>                | 0.05% Sodium Azide  |
| <b>Isotype</b>                     | IgG   |
| <b>Conjugate</b>                   | DyLight 650   |
| <b>Purity</b>                      | Immunogen affinity purified   |
| <b>Buffer</b>                      | 50mM Sodium Borate  |
| <b>Product Description</b>         |   |
| <b>Host</b>                        | Rabbit  |
| <b>Gene ID</b>                     | 64344   |
| <b>Gene Symbol</b>                 | HIF3A   |
| <b>Species</b>                     | Human, Mouse  |
| <b>Immunogen</b>                   | Bacterially expressed mouse HIF-3 alpha (N-terminus). [UniProt# Q0VBL6]   |
| <b>Notes</b>                       | Dylight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.   |
| <b>Product Application Details</b> |   |
| <b>Applications</b>                | Western Blot  |
| <b>Recommended Dilutions</b>       | Western Blot  |
| <b>Application Notes</b>           | <p>This HIF-3 alpha antibody is useful for Western blot, where a band is observed ~70 kDa.</p> <p>The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.</p> |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-2287C](http://www.novusbio.com/reviews/submit/NB100-2287C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

