

# Product Datasheet

## MINDY3 Antibody - BSA Free

### NB100-2504

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

#### Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-2504](http://www.novusbio.com/NB100-2504)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-2504](http://www.novusbio.com/reviews/destination/NB100-2504)



**NB100-2504**

MINDY3 Antibody - BSA Free

**Product Information**

|                      |                                                                                        |
|----------------------|----------------------------------------------------------------------------------------|
| <b>Unit Size</b>     | 0.1 ml                                                                                 |
| <b>Concentration</b> | 1 mg/ml                                                                                |
| <b>Storage</b>       | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| <b>Clonality</b>     | Polyclonal                                                                             |
| <b>Preservative</b>  | 0.02% Sodium Azide                                                                     |
| <b>Isotype</b>       | IgG                                                                                    |
| <b>Purity</b>        | Immunogen affinity purified                                                            |
| <b>Buffer</b>        | PBS                                                                                    |

**Product Description**

|                         |                                                                                                                                                                                                                                         |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>      | Novus Biologicals Rabbit MINDY3 Antibody - BSA Free (NB100-2504) is a polyclonal antibody validated for use in ELISA. Anti-MINDY3 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| <b>Host</b>             | Rabbit                                                                                                                                                                                                                                  |
| <b>Gene ID</b>          | 80013                                                                                                                                                                                                                                   |
| <b>Gene Symbol</b>      | MINDY3                                                                                                                                                                                                                                  |
| <b>Species</b>          | Human                                                                                                                                                                                                                                   |
| <b>Reactivity Notes</b> | Human                                                                                                                                                                                                                                   |
| <b>Immunogen</b>        | Residues 59-76 AFLKKLLFSSEKSSWRD of human MINDY3 protein.                                                                                                                                                                               |

**Product Application Details**

|                              |                                                                                                                                                                                                                                                                   |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>          | ELISA                                                                                                                                                                                                                                                             |
| <b>Recommended Dilutions</b> | ELISA 1:100-1:2000                                                                                                                                                                                                                                                |
| <b>Application Notes</b>     | This antibody is a high-titer, purified antibody that has been tested only by ELISA. while it is expected that this antibody will preform well in Western Blot and IHC applications, there is no guarantee that it will work in any application other than ELISA. |

**Publications**

Liu B, Liu Y, Chen J et al. CARP is a novel caspase recruitment domain containing pro-apoptotic protein. Biochem Biophys Res Commun. 2002-01-01 [PMID: 12054670]





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-2504**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00080013-P01-10ug | Recombinant Human MINDY3 GST (N-Term) Protein       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-2504](http://www.novusbio.com/reviews/submit/NB100-2504)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

