

Product Datasheet

HIF-3 alpha Antibody [DyLight 550] NB100-2529R

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-2529R

Updated 9/10/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-2529R



NB100-2529R

HIF-3 alpha Antibody [DyLight 550]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 550
Purity	IgG purified
Buffer	50mM Sodium Borate
Target Molecular Weight	68 kDa
Product Description	
Host	Rabbit
Gene ID	64344
Gene Symbol	HIF3A
Species	Human, Mouse
Immunogen	A synthetic peptide made to an internal region of human HIF-3 alpha (within residues 250-350). [Swiss-Prot: Q66K72]
Notes	Dylight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation
Application Notes	<p>This HIF-3 alpha antibody is useful for Immunoprecipitation, Western Blot and Immunocytochemistry/Immunofluorescence. In Western Blot analysis, a band is recognized at ~68 kDa representing HIF-3 alpha. Because a band at ~71 kDa is also seen in hypoxic samples, it is important to run a negative control lysate as well (see image). In ICC/IF, cytoplasmic and nuclear staining was observed in HeLa cells, which is consistent with the literature under normoxic conditions for HIF-3 alpha (PMID: 16775626).</p> <p>The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.</p>



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-2529R

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

