

Product Datasheet

FANCC Antibody NB100-2565

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-2565

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-2565



NB100-2565**FANCC Antibody****Product Information**

Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Unpurified
Buffer	Whole antisera
Target Molecular Weight	65 kDa

Product Description


Host	Rabbit
Gene ID	2176
Gene Symbol	FANCC
Species	Human
Reactivity Notes	Immunogen displays the following percentage of sequence identity for non-tested species: rat (92%), mouse (92%) and bovine (81%).
Immunogen	A C-terminal peptide within residues 450-558 of human FANCC. [SwissProt# Q00597]

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1:1000
Application Notes	This FANCC antibody is useful for Western Blot where a band is seen at approximately 65 kDa. The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.

Images

Western Blot: FANCC Antibody [NB100-2565] - Detection of FANCC in transfected COS1 cell lysate.



Publications

Yung SK, Tilgner K, Ledran MH et al. Human pluripotent stem cell models of Fanconi Anaemia deficiency reveal an important role for Fanconi Anaemia proteins in cellular reprogramming and survival of haematopoietic progenitors. *Stem Cells* 2012-12-27 [PMID: 23280624] (WB)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-2565

NBL1-10587	FANCC Overexpression Lysate
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-2565

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

