Product Datasheet

SRC1 Antibody - BSA Free NB100-312

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com

technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NB100-312

Updated 2/21/2025 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NB100-312



NB100-312

SRC1 Antibody - BSA Free

Product Information	
Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)
Product Description	
Description	This antibody can be used as the primary antibody in a PLA assay with the following as complementing antibodies:NB500-160, NBP1-19102
Host	Rabbit
Gene ID	8648
Gene Symbol	NCOA1
Species	Human, Mouse
Immunogen	A 29 amino acid synthetic peptide, which represented a portion of human nuclear receptor coactivator 1 encoded within exon 12 (LocusLink ID 8648).
Notes	This antibody can be used as the primary antibody in a PLA assay with the following as complementing antibodies:NB500-160, NBP1-19102
Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunohistochemistry, Immunoprecipitation
Application Notes	NB 100-312 may be used for Western blot. The investigator should determine the optimal working dilution for a specific application.







www.novusbio.com



technical@novusbio.com

Serum protocol for SRC1 Antibody (NB100-312)

Nuclear Extract and Cytoplasmic Fraction Preparation protocol for SRC1 Antibody (NB100-312): Nuclear Extract and Cytoplasmic Fraction Preparation

1. Nuclear extracts (NE) and cytoplasmic fractions (S100) were prepared by Dignam's method (Dignam, Lebovitz, and Roeder, Nucleic Acids Res. 11: 1475-1489. 1983).

2. 100 liters of HeLa cell culture were harvested and washed 3 times with cold PBS.

3. The packed-cell volume (PCV) was measured, and the cell pellet was gently resuspended with 5 PCVs of hypotonic buffer (10 mM HEPES-KOH [pH 8], 10 mM KCl, 1.5 mM MgCl2, 1 mM DTT, 0.2 mM PMSF).

4. Cells were incubated on ice for 10 minutes and then pelleted by centrifugation at 1,800xg for 10 minutes.

5. Hypotonic buffer was added to 2 PCVs, and cells were resuspended and then homogenized with 15 strokes using a pestle B in a Dounce glass homogenizer until the cells were more than 90% lysed, as determined by a light microscope.

6. The lysate was centrifuged at 20,000xg for 30 minutes at 4 degrees Celcius.

7. The supernatant was saved for S100 fraction, and the pellet was saved to measure the packed nuclear volume (PNV).

8. 0.4 ml of extraction buffer (20 mM HEPES-KOH [pH 8], 0.6 M KCl, 1.5 mM MgCl2, 0.2 mM EDTA, 25% [vol/vol] glycerol, 1 mM DTT, 0.2 mM PMSF) per ml of PNV was added.

9. Cell nuclei were homogenized with 10 strokes of pestle A in the homogenizer.

10. Suspension was stirred at 4 degrees Celcius for 30 minutes and centrifuged for 30 minutes at 20,000xg.

11. The supernatant (nuclear extract) was aliquotted for use.

12. The S100 fraction (resulting supernatant) was mixed with 0.11 volume of high-salt buffer (20 mM HEPES-KOH [pH 8], 1.2 M KCl, 1.5 mM MgCl2, 0.2 mM EDTA, 20% [vol/vol] glycerol, 1 mM DTT, 0.2 mM PMSF) and centrifuged at 100,000xg for 60 minutes at 4 degrees Celcius.

13. This supernatant was dialyzed for 2 hours at 4 degrees Celcius.

14. The sample was centrifuged for 30 minutes at 20,000xg and the supernatant (S100) was aliquotted for use.

Immunoprecipitation

Antibody characterization:

1. HeLa NE and S100 were diluted with 1 volume of RIPA buffer [150 mM NaCl, 1% NP-40, 0.5% DOC, 0.1% SDS, 50 mM Tris [pH 8]).

2. Cleared by spinning at 100,000 g for 20 minutes at 4 degrees Celcius.

3. 1 ml of supernatant (~10 mg total protein) was mixed with 20 ug of primary antibody (NB 100-312) and rotated overnight at 4 degrees Celcius.

4. Supernatant was mixed with 0.05 ml of protein A-sepharose beads (50% slurry) and rotated for 2 hours at 4 degrees Celcius.

5. Immunoprecipitates were washed 3 times with the 10% RIPA in PBS.

6. The washed beads were boiled with 0.04 ml of Laemmli buffer and subjected to SDS-PAGE (4-20% Tris-glycine gel).

Complex purification:

1. NE and S100 were cleared by spinning at 20,000 g for 30 minutes at 4 degrees Celcius.

2. 1.5 ml of supernatant (~15 mg total protein) was mixed with 20 ug of primary antibody (NB 100-312) and rotated for 4 hours at 4 degrees Celcius.

3. Sample and antibody mixture were centrifuged at 15,000 g for 20 minutes at 4 degrees Celcius.

4. Supernatant was mixed with 0.05 ml of protein A-sepharose beads (50% slurry) and rotated for 1 hour at 4 degrees Celcius.

5. Immunoprecipitates were washed 3 times with the NETN buffer (20 mM Tris-HCI [pH 8], 100 mM NaCI, 1 mM EDTA, 0.5% NP-40).

6. The washed beads were boiled with 0.04 ml of Laemmli buffer and subjected to SDS-PAGE (4-20% Tris-glycine gel).

*If an insufficient amount of protein is purified for identification from 15 mg of extract, carry out the same procedure using 50-100 mg of extract to increase the amount of purified protein yield.

www.novusbio.com





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

Products Related to NB100-312

NB800-PC9	HeLa Nuclear Cell Lysate
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-312

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

