

# Product Datasheet

## SAPS1 Antibody NB100-41413

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB100-41413](http://www.novusbio.com/NB100-41413)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB100-41413](http://www.novusbio.com/reviews/destination/NB100-41413)



**NB100-41413****SAPS1 Antibody**


| Product Information |                             |
|---------------------|-----------------------------|
| Unit Size           | 0.1 ml                      |
| Concentration       | 0.2 mg/ml                   |
| Storage             | Store at 4C. Do not freeze. |
| Clonality           | Polyclonal                  |
| Preservative        | 0.09% Sodium Azide          |
| Isotype             | IgG                         |
| Purity              | Immunogen affinity purified |
| Buffer              | TBS and 0.1% BSA            |

| Product Description |                                                                                                                                                                                                                      |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Host                | Rabbit                                                                                                                                                                                                               |
| Gene ID             | 22870                                                                                                                                                                                                                |
| Gene Symbol         | PPP6R1                                                                                                                                                                                                               |
| Species             | Human                                                                                                                                                                                                                |
| Immunogen           | The immunogen recognized by this antibody maps to a region between residue 841 and the C-terminus (residue 891) of human SAPS domain family, member 1 using the numbering given in entry NP_055746.2 (GenelD 22870). |

| Product Application Details |                                                                   |
|-----------------------------|-------------------------------------------------------------------|
| Applications                | Western Blot, Immunoprecipitation                                 |
| Recommended Dilutions       | Western Blot 1:2000-1:10000, Immunoprecipitation 2-5 ug/mg lysate |

**Images**

Western Blot: SAPS1 Antibody [NB100-41413] - Detection of Human SAPS1 on HeLa whole cell lysate using NB100-41413. SAPS1 was also immunoprecipitated by another rabbit anti-SAPS1 antibody.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB100-41413**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |
| NBP2-56076PEP | SAPS1 Recombinant Protein Antigen                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB100-41413](http://www.novusbio.com/reviews/submit/NB100-41413)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

