

Product Datasheet

SCAF11/SFRS2IP Antibody

NB100-55266

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Reviews: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-55266

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-55266



NB100-55266**SCAF11/SFRS2IP Antibody**

Product Information	
Unit Size	100 ul
Concentration	0.2 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	TBS and 0.1% BSA


Product Description	
Description	Not to be confused with rodent Caspase 11, which is synonymouse with Caspase 4. Human Caspase 11 alternate names: anti-arginine/serine-rich 2 antibody, anti-CASP11 antibody, anti-Caspase 11 apoptosis related cysteine protease, anti-caspase 4 antibody, anti-Caspase 4 apoptosis related cysteine peptidas, anti-Caspase11 antibody, Anti-Caspase4 antibody, anti-CASPL antibody, anti-CTD-associated SR protein 11 antibody, anti-ICH 3 antibody, anti-ICH3 antibody, anti-ICH3 protease antibody, anti-interacting protein antibody, anti-SC35-interacting protein 1 antibody, anti-SFRS2IP antibody, anti-SIP1 antibody, anti-splicing factor antibody, anti-Splicing regulatory protein 129 antibody, anti-SRrp129 antibody
Host	Rabbit
Gene ID	9169
Gene Symbol	SCAF11
Species	Human
Immunogen	The immunogen recognized by this antibody maps to a region between residue 1098 and the C-terminus (residue 1148) of human splicing factor, arginine/serine-rich 2, interacting protein using the numbering given in entry CAA72121.1 (GeneID 9169).
Notes	Not to be confused with rodent Caspase 11, which is synonymouse with Caspase 4. Human Caspase 11 alternate names: anti-arginine/serine-rich 2 antibody, anti-CASP11 antibody, anti-Caspase 11 apoptosis related cysteine protease, anti-caspase 4 antibody, anti-Caspase 4 apoptosis related cysteine peptidas, anti-Caspase11 antibody, Anti-Caspase4 antibody, anti-CASPL antibody, anti-CTD-associated SR protein 11 antibody, anti-ICH 3 antibody, anti-ICH3 antibody, anti-ICH3 protease antibody, anti-interacting protein antibody, anti-SC35-interacting protein 1 antibody, anti-SFRS2IP antibody, anti-SIP1 antibody, anti-splicing factor antibody, anti-Splicing regulatory protein 129 antibody, anti-SRrp129 antibody

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000-1:10000, Immunohistochemistry 1:200 - 1:1000, Immunocytochemistry/ Immunofluorescence 1:100 to 1:500, Immunoprecipitation 2-5 ug/mg lysate, Immunohistochemistry-Paraffin 1:200 - 1:1000
Application Notes	Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections. SCAF11/SFRS2IP atibody validated for ICC/IF from a customer review.




Images

Western Blot: SCAF11/SFRS2IP Antibody [NB100-55266] - Detection of Human SFRS2IP on HeLa whole cell lysate using NB100-55266. SFRS2IP was also immunoprecipitated by another rabbit anti-SFRS2IP antibody.



Immunohistochemistry-Paraffin: SCAF11/SFRS2IP Antibody [NB100-55266] - Sample: FFPE section of human metastatic lymph node. Antibody: Affinity purified rabbit anti- SFRS2IP used at a dilution of 1:200 (1ug/ml). Detection: DAB





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-55266

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-55306PEP	SCAF11/SFRS2IP Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-55266

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

