

Product Datasheet

BAG2 Antibody NB100-56087

Unit Size: 0.05 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-56087

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-56087




NB100-56087**BAG2 Antibody**

Product Information	
Unit Size	0.05 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Unpurified
Buffer	Whole antisera
Product Description	
Host	Rabbit
Gene ID	9532
Gene Symbol	BAG2
Species	Human
Immunogen	Full-length recombinant human BAG-2 protein was used as immunogen.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation
Recommended Dilutions	Western Blot 1:1000-1:2000, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:1000-1:5000, Immunoprecipitation 1:50-1:200, Immunohistochemistry-Paraffin 1:1000-1:5000




Images


Western Blot: BAG2 Antibody [NB100-56087] - analysis of BAG-2 using this antibody. HeLa stably transfected Myc-tagged CHIP (



Immunocytochemistry/Immunofluorescence: BAG2 Antibody [NB100-56087] - Immunofluorescence confocal microscopy of BAG-2 using NB100-56087 at 1:2000. HeLa cells stably transfected with Myc-CHIP (carboxyl terminus of Hsp70-interacting protein) and growing on coverslips were left untreated (A) or heat shocked at 42 degrees for 30 min (B). Cells were fixed with 3.7% paraformaldehyde, permeabilized with 0.5% Triton X-100 prior to antibody staining. BAG-2 localized to the cytoplasm of untreated cells, and to the nucleus of cells subjected to heat shock. Please refer to Dai et al (2005) for additional information on HeLa cells stably transfected with Myc-CHIP.



Immunohistochemistry-Paraffin: BAG2 Antibody [NB100-56087] - BAG-2 was detected in immersion fixed paraffin-embedded sections of human testis using Rabbit Anti-Human BAG-2 polyclonal Antibody (Catalog # NB100-56087) at 1:1000 for 1 hour at room temperature followed by incubation with the Anti-Rabbit IgG VisUCyte™ HRP Polymer Antibody (Catalog # VC003). Tissue was stained using DAB (brown) and counterstained with hematoxylin (blue). Specific staining was localized to the cytoplasm and plasma membrane.



Publications

Dai Q, Qian SB, Li HH et al. Regulation of the cytoplasmic quality control protein degradation pathway by BAG2. J Biol Chem 2005-11-18 [PMID: 16169850] (ICC/IF, WB)

Wu F, Wang P, Young LC et al. Proteome-Wide Identification of Novel Binding Partners to the Oncogenic Fusion Gene Protein, NPM-ALK, using Tandem Affinity Purification and Mass Spectrometry. Am J Pathol 2009-01-08 [PMID: 19131589]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-56087

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-51025-0.1mg	Recombinant Human BAG2 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-56087

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

