

Product Datasheet

IFN-alpha 2 Antibody - BSA Free NB100-75556

Unit Size: 0.05 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-75556

Updated 7/16/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-75556



NB100-75556

IFN-alpha 2 Antibody - BSA Free

Product Information

Unit Size	0.05 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2)

Product Description

Host	Chicken
Gene ID	3440
Gene Symbol	IFNA2
Species	Human
Reactivity Notes	This antibody is reactive against Human.
Specificity/Sensitivity	Interferon alpha 2b
Immunogen	The immunogen for this antibody was a sequence from 25a.a. to 186a.a. from Interferon alpha 2b(a 188a.a. long protein)
Notes	Chicken products cannot be exported to Canada.


Product Application Details

Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 1:200-1:500, ELISA 1:200-1:500
Application Notes	ELISA: Recombinant protein as test antigen. Affinity-Purified anti- IFNa2b IgY as primary antibody (0.1 ug/ml). and Goat anti-IgY-HRP as 2nd antibody.




Images

Western Blot: Recombinant human IFNa2b as antigen. IgY dilution: 1:2000; Goat anti-IgY Fc-HRP: 1: 1000. Colorimetric method for signal development. (anti-IFNA2)



ELISA: Interferon alpha 2b Antibody [NB100-75556] - Recombinant protein as test antigen. Affinity-Purified anti- IFNa2b IgY as primary antibody (0.1 ug/ml), and Goat anti-IgY-HRP as 2nd antibody. (anti-IFNA2)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-75556

NBL1-11842	IFN-alpha 2 Overexpression Lysate
BAF010	Goat anti-Chicken IgY Secondary Antibody [Biotin]
NB7276	Goat anti-Chicken IgM Heavy Chain Secondary Antibody
NBP2-34971-100ug	Recombinant Human IFN-alpha 2 (alpha 2a) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-75556

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

