

Product Datasheet

STAMP2/STEAP4 Antibody - BSA Free NB100-80831

Unit Size: 0.05 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-80831

Updated 4/3/2022 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-80831



NB100-80831

STAMP2/STEAP4 Antibody - BSA Free

Product Information

Unit Size	0.05 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS

Product Description

Host	Rabbit
Gene ID	79689
Gene Symbol	STEAP4
Species	Human, Mouse, Rat
Specificity/Sensitivity	At least two isoforms are known to exist for this protein. This STAMP2/STEAP4 antibody does not cross-react with other STEAP proteins.
Immunogen	A 14 aa synthetic peptide near the center of human STAMP2/STEAP4. The immunogen is located within aa 210 - 260 of STAMP2/STEAP4.


Product Application Details

Applications	Western Blot, ELISA, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 2 ug/ml, ELISA, Immunohistochemistry 5 ug/ml, Immunohistochemistry-Paraffin 5 ug/ml
Application Notes	WB: Successful result is reported by anonymous customer.




Images

Western Blot: STAMP2/STEAP4 Antibody [NB100-80831] - Rat liver tissue lysate.



Immunohistochemistry-Paraffin: STAMP2/STEAP4 Antibody [NB100-80831] - STEAP4 Antibody [NB100-80831] - Staining of human heart formalin-fixed, paraffin-embedded tissue after heat-induced antigen retrieval. Antibody concentration 5 ug/ml.



Publications

Chuang CT, Guh JY, Lu CY et al. Steap4 attenuates high glucose and S100B-induced effects in mesangial cells. J. Cell. Mol. Med. 2015 Mar 27 [PMID: 25817898] (IHC-P, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NB100-80831

NBL1-16537	STAMP2/STEAP4 Overexpression Lysate
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-80831

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

