

# Product Datasheet

## Adiponectin/Acrp30 Antibody NB110-6026

Unit Size: 0.2 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB110-6026](http://www.novusbio.com/NB110-6026)

Updated 5/28/2019 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB110-6026](http://www.novusbio.com/reviews/destination/NB110-6026)



**NB110-6026****Adiponectin/Acrp30 Antibody**

<b>Product Information</b>	
<b>Unit Size</b>	0.2 mg
<b>Concentration</b>	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.05% Sodium Azide
<b>Purity</b>	Immunogen affinity purified
<b>Buffer</b>	PBS, 1 mg/ml BSA

<b>Product Description</b>	
<b>Host</b>	Rabbit
<b>Gene ID</b>	9370
<b>Gene Symbol</b>	ADIPOQ
<b>Species</b>	Mouse, Human (Negative)
<b>Reactivity Notes</b>	Cross-reacts with Mouse and very weakly with Rat. Does not cross-react with Human. Not yet tested in other species. Please see the enquiries section of the online datasheet for further discussion.
<b>Specificity/Sensitivity</b>	Detects Acrp30 from mouse serum.
<b>Immunogen</b>	Synthetic Peptides: E(18)DDVTTEELAPALV(32) and F(187)TYDQYQEKNVDQA(200). Immunizing peptides correspond to amino acid residues 18-32 and 187-200 from mouse Acrp 30 protein.

<b>Product Application Details</b>	
<b>Applications</b>	Western Blot, Immunoprecipitation
<b>Recommended Dilutions</b>	Western Blot, Immunoprecipitation
<b>Application Notes</b>	Western Blot: use at a concentration of 4-8 ug/ml. Immunoprecipitation: assay dependent. By Western blot, this antibody detects an ~30 kDa protein representing Acrp 30 from mouse serum.

**Images**

Western Blot: Adiponectin Antibody [NB110-6026] - Detection of Adiponectin from mouse serum.

**Fig. 1**



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Novus Biologicals Canada**

461 North Service Road West, Unit B37  
Oakville, ON L6M 2V5  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada@novusbio.com

### **Novus Biologicals Europe**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB110-6026](http://www.novusbio.com/reviews/submit/NB110-6026)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

