

Product Datasheet

PDGF R beta Antibody (MM0005-5C37) NB110-60970

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-60970

Updated 7/15/2019 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-60970



NB110-60970

PDGF R beta Antibody (MM0005-5C37)

Product Information	
Unit Size	0.1 mg
Concentration	LYOPH mg/ml
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	MM0005-5C37
Preservative	No Preservative
Reconstitution Instructions	Reconstitute with 0.5 ml sterilized PBS to a final concentration of 0.2 mg/ml.
Isotype	IgG1
Purity	Protein G purified
Buffer	A 0.2 um filtered solution in PBS


Product Description	
Host	Mouse
Gene ID	5159
Gene Symbol	PDGFRB
Species	Human, Rabbit
Reactivity Notes	Rabbit reactivity reported in scientific literature (PMID: 30125282).
Specificity/Sensitivity	This antibody detects human PDGFRb in western blots. No cross-reactivity is shown with hPDGFRa
Immunogen	Recombinant hPGFDR beta N-terminal

Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, CyTOF-ready
Recommended Dilutions	Western Blot 1:400-1:1000, Flow Cytometry 12-5 ul/10 6 cells , Immunohistochemistry 1:50 -500, Immunocytochemistry/Immunofluorescence, Immunohistochemistry-Paraffin 1:50-500, CyTOF-ready
Application Notes	This antibody is Cytof ready. Use in Immunocytochemistry/immunofluorescence reported in scientific literature (PMID: 30125282).




Images

Immunohistochemistry-Paraffin: PDGF Receptor beta Antibody (MM0005-5C37) [NB110-60970] - Formalin fixed and paraffin embedded normal human placenta tissue is subjected to mouse anti human PDGFR beta



Immunohistochemistry-Paraffin: PDGF Receptor beta Antibody (MM0005-5C37) [NB110-60970] - Formalin fixed and paraffin embedded normal human placenta tissue is subjected to mouse anti human PDGFR beta.



Publications

Szubert, S;Moszynski, R;Szpurek, D;Romaniuk, B;Sajdak, S;Nowicki, M;Michalak, S; The expression of Platelet-Derived Growth Factor Receptors (PDGFRs) and their correlation with overall survival of patients with ovarian cancer Ginekol. Pol. Apr 12 2019 12:00AM [PMID: 30977516] (IHC, Human)

Mao Y, Liu X, Song Y et al. VEGF-A/VEGFR-2 and FGF-2/FGFR-1 but not PDGF-BB/PDGFR-b play important roles in promoting immature and inflammatory intraplaque angiogenesis PLoS ONE Aug 20 2018 12:00AM [PMID: 30125282] (IHC, ICC/IF, WB, Rabbit)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-60970

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

