

Product Datasheet

TIMP-4 Antibody (MM0037-9F3) - Azide and BSA Free NB110-61003

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-61003

Updated 4/7/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-61003



NB110-61003

TIMP-4 Antibody (MM0037-9F3) - Azide and BSA Free

Product Information

| | |
|------------------------------------|---|
| Unit Size | 0.1 mg |
| Concentration | LYOPH mg/ml |
| Storage | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | MM0037-9F3 |
| Preservative | No Preservative |
| Reconstitution Instructions | Reconstitute with sterilized PBS to a final concentration of 0.2 mg/ml. |
| Isotype | IgG1 |
| Purity | Protein G purified |
| Buffer | Lyophilized from a 0.2 um filtered solution in PBS. 0.025 mg size is provided in liquid form, PBS |

Product Description

| | |
|--------------------------------|---|
| Host | Mouse |
| Gene ID | 7079 |
| Gene Symbol | TIMP4 |
| Species | Human |
| Specificity/Sensitivity | No cross-reactivity to human TIMP-1, 2 and 3 is noted |
| Immunogen | Human recombinant TIMP-4 |

Product Application Details

| | |
|------------------------------|--|
| Applications | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Western Blot 1:100-1:1000, Immunohistochemistry 1:100-200, Immunohistochemistry-Paraffin 1:100-1:200 |

Images

Immunohistochemistry: TIMP-4 Antibody (MM0037-9F3) [NB110-61003]
 - Formalin fixed and paraffin embedded normal human placenta tissue
 using NB110-61003





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-61003

| | |
|------------------|---|
| NBL1-16927 | TIMP-4 Overexpression Lysate |
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB720-B | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-97005-0.5mg | Mouse IgG1 Isotype Control (MG1) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-61003

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

