

Product Datasheet

ARAP1 Antibody (ARAP1-2) - BSA Free NB110-68801

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-68801

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-68801



NB110-68801

ARAP1 Antibody (ARAP1-2) - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	ARAP1-2
Preservative	0.9% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.4)
Target Molecular Weight	160 kDa

Product Description


Host	Mouse
Gene ID	116985
Gene Symbol	ARAP1
Species	Human
Specificity/Sensitivity	The ARAP1-2 reacts with C-terminal part of ARAP1, a 160 kDa adaptor protein.
Immunogen	His6-tagged protein corresponding to amino acids 1190-1450 of human ARAP1

Product Application Details

Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 2 ug/ml, Immunocytochemistry/ Immunofluorescence 1:10-1:2000, Immunoprecipitation 1:50

Images

Immunocytochemistry/Immunofluorescence: ARAP1 Antibody (ARAP1-2) [NB110-68801] - Colocalization of ARAP1 with DR4 at the plasma membrane and in early endosomes. ARAP1-transfected NCTC cells were treated with TRAIL for 0, 30 or 60 and fixed with methanol. DR4 was detected by rabbit polyclonal (red) and ARAP1 with ARAP1-2 mouse monoclonal antibody (green). Arrowheads show ARAP1-DR4 colocalization at the plasma membrane and intracellular membranes, whereas arrows show their colocalization in early endosomes.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-68801

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NB110-68801UV	ARAP1 Antibody (ARAP1-2) [DyLight 350]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-68801

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

