

Product Datasheet

PRMT3 Antibody (PRMT3-367) - BSA Free NB120-12190

Unit Size: 0.05 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB120-12190

Updated 6/13/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB120-12190



NB120-12190


PRMT3 Antibody (PRMT3-367) - BSA Free

Product Information	
Unit Size	0.05 ml
Concentration	2.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	PRMT3-367
Preservative	0.9% Sodium Azide
Isotype	IgG3
Purity	Protein A or G purified
Buffer	10mM PBS (pH 7.4)
Product Description	
Host	Mouse
Gene ID	10196
Gene Symbol	PRMT3
Species	Human, Mouse, Rat
Reactivity Notes	Not yet tested in other species. Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information.
Specificity/Sensitivity	does not cross react with PRMT1, 2, 4, 5, or 6.
Immunogen	Synthetic peptide: TVHKNKKDPRSLIV conjugated to KLH, corresponding to amino acids 501-514 of Rat PRMT3.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:100-1:2000, ELISA 1:100-1:2000, Immunocytochemistry/ Immunofluorescence 1:10-1:2000
Application Notes	ELISA: Use at an assay dependent dilution. ICC: Use at an assay dependent dilution. Western Blot: Use at a concentration of 1-2 ug/ml. This concentration is recommended using a whole cell extract of the 293T cell line transfected with human PRMT3. Predicted molecular weight: 58 kDa. Optimal dilutions/concentrations should be determined by the end user.



Images

Western Blot: PRMT3 Antibody (PRMT3-367) - BSA Free [NB120-12190] - Whole cells extract of 293T cell line transfected with human PRMT3 were separated on SDS-PAGE and probed with Monoclonal Anti-PRMT3 antibody produced in Mouse, Clone: PRMT3-367. The antibody was developed using 1:5000 Anti-Mouse IgG (Fab specific)- Peroxidase antibody produced in Goat.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB120-12190

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

