

# Product Datasheet

## S5a/Angiocidin Antibody NB120-3318

Unit Size: 200uL

Store at -20C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB120-3318](http://www.novusbio.com/NB120-3318)

Updated 5/15/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB120-3318](http://www.novusbio.com/reviews/destination/NB120-3318)



**NB120-3318****S5a/Angiocidin Antibody**


| <b>Product Information</b> |  |
|----------------------------|--|
| <b>Unit Size</b>           | 200uL                                    |
| <b>Concentration</b>       | 1 mg/ml                                  |
| <b>Storage</b>             | Store at -20C. Avoid freeze-thaw cycles. |
| <b>Clonality</b>           | Polyclonal                               |
| <b>Preservative</b>        | 0.05% Sodium Azide                       |
| <b>Isotype</b>             | IgG                                      |
| <b>Purity</b>              | Ammonium sulfate precipitation           |
| <b>Buffer</b>              | PBS with 1 mg/ml BSA                     |

| <b>Product Description</b>     |   |
|--------------------------------|---|
| <b>Host</b>                    | Rabbit  |
| <b>Gene ID</b>                 | 5710  |
| <b>Gene Symbol</b>             | PSMD4   |
| <b>Species</b>                 | Human   |
| <b>Specificity/Sensitivity</b> | Detects proteasome 19S subunit S5A.           |
| <b>Immunogen</b>               | Human recombinant proteasome 19S subunit S5A. |

| <b>Product Application Details</b> |  |
|------------------------------------|--|
| <b>Applications</b>                | Western Blot, Immunoprecipitation  |
| <b>Recommended Dilutions</b>       | Western Blot 1:500, Immunoprecipitation 1:10 - 1:500   |
| <b>Application Notes</b>           | WB: Detects an approx. 50 kDa protein representing proteasome 19S subunit S5A in HeLa cell lysate. |

**Images**

Western Blot: Proteasome 19S S5A Antibody [NB120-3318] - Analysis of HeLa cell lysate. **Fig. 1**





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB120-3318](http://www.novusbio.com/reviews/submit/NB120-3318)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

