

Product Datasheet

p19ARF/CDKN2A Antibody (5-C3-1) [FITC] NB200-174F

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB200-174F

Updated 10/23/2024 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB200-174F



NB200-174F

p19ARF/CDKN2A Antibody (5-C3-1) [FITC]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	5-C3-1
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	FITC
Purity	Protein A purified
Buffer	PBS
Product Description	
Host	Rat
Gene ID	12578
Gene Symbol	Cdkn2a
Species	Human, Mouse, Golden Syrian Hamster (Negative)
Reactivity Notes	Human reactivity reported in scientific literature (PMID: 27177758). The antibody does not react with Golden Syrian hamster. Armenian hamster has not been tested.
Immunogen	A synthetic peptide containing amino acids 54-75 of the murine p19ARF/CDKN2A. [UniProt# Q64364]
Product Application Details	
Applications	Western Blot, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, Immunoprecipitation, CyTOF-ready, Knockout Validated
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Flow (Intracellular), CyTOF-ready, Knockout Validated
Application Notes	Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB200-174F

NB200-174H	p19ARF/CDKN2A Antibody (5-C3-1) [HRP]
NB200-106PEP	p19ARF/CDKN2A Antibody Blocking Peptide
NBP3-38956	Mouse p19ARF/CDKN2A ELISA Kit (Colorimetric)
NB200-103	p53 Antibody (PAb 240) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB200-174F

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

