

Product Datasheet

p19ARF/CDKN2A Antibody (5-C3-1) [7C] NB200-174V2

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB200-174V2

Updated 7/28/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB200-174V2




NB200-174V2

p19ARF/CDKN2A Antibody (5-C3-1) [7C]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	5-C3-1
Preservative	0.05% Sodium Azide
Isotype	IgG2b
Conjugate	7C
Purity	Protein A purified
Buffer	PBS
Target Molecular Weight	19 kDa
Product Description	
Host	Rat
Gene ID	12578
Gene Symbol	Cdkn2a
Species	Human, Mouse, Golden Syrian Hamster (Negative)
Reactivity Notes	Human reactivity reported in scientific literature (PMID: 27177758). The antibody does not react with Golden Syrian hamster. Armenian hamster has not been tested.
Immunogen	A synthetic peptide containing amino acids 54-75 of the murine p19ARF/CDKN2A. [UniProt# Q64364]
Notes	This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen
Application Notes	<p>This p19ARF/CDKN2A (5-C3-1) antibody is useful for Immunocytochemistry/Immunofluorescence, Immunohistochemistry on frozen sections, Immunoprecipitation, Flow Cytometry and Western blot, where a band can be seen at approximately 19 kDa.</p> <p>The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.</p>

Images

Flow Cytometry: p19ARF/CDKN2A Antibody (5-C3-1) [7C] - Conducted the antibody trial with two different dilutions of p19ARF/CDKN2A Antibody (5-C3-1) [PE] (Catalog # NB200-174PE) (1:500 and 1:100) and compared them to an IgG control which was also PE labelled. Image from a verified customer review.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB200-174V2

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

