

Product Datasheet

Annexin A10 Antibody NB300-894

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-894

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-894



NB300-894**Annexin A10 Antibody****Product Information**

Unit Size	0.1 mg
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA

Product Description


Host	Goat
Gene ID	11199
Gene Symbol	ANXA10
Species	Human
Immunogen	Peptide with sequence FCGDYVQGTIFPAPN corresponding to N-Terminus according to NP_009124.2.

Product Application Details

Applications	Immunohistochemistry, Immunohistochemistry-Paraffin, Peptide ELISA
Recommended Dilutions	Immunohistochemistry 3 - 5 ug/ml, Immunohistochemistry-Paraffin 3 - 5 ug/ml, Peptide ELISA Detection limit 1:8000
Application Notes	WB: Preliminary experiments (at 1 ug/ml) gave an approx. 45-48 kDa band in Hela, A431 and human kidney lysates. Please note that currently we cannot find an explanation in scientific literature for the band we observe given the predicted MW of approx. 37 kDa band according to NP_009124. IHC-P: Human pancreas shows staining of islet of Langerhans.

Images

Immunohistochemistry-Paraffin: Annexin A10 Antibody [NB300-894] - Human Pancreas using this antibody at 3.8ug/ml. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.

**Publications**

Morgan RO, Jenkins NA, Gilbert DJ et al. Novel human and mouse annexin A10 are linked to the genome duplications during early chordate evolution. Genomics 1999-08-15 [PMID: 10458909]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-894

NBL1-07558	Annexin A10 Overexpression Lysate
HAF017	Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
HAF109	Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB410-28088-1mg	Goat IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-894

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



