Product Datasheet

ABCA1 Antibody (3A1.891.3) NB400-164

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 9

Protocols, Publications, Related Products, Reviews, Research Tools and Images at: www.novusbio.com/NB400-164

Updated 10/23/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications Submit a review at www.novusbio.com/reviews/destination/NB400-164



NB400-164

ABCA1 Antibody (3A1.891.3)

Product Information					
Unit Size	0.1 ml				
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.				
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.				
Clonality	Monoclonal				
Clone	3A1.891.3				
Preservative	0.1% Sodium Azide				
Isotype	IgG1 Kappa				
Purity	Unpurified				
Buffer	Ascites				
Target Molecular Weight	250 kDa				
Product Description					
Host	Rat				
Gene ID	19				
Gene Symbol	ABCA1				
Species	Mouse, Human (Negative)				
Immunogen	A1-N3 peptide (within residues 200-250) of mouse ABCA1 Antibody (3A1.891.3) covalently linked to chicken egg albumin by means of glutaraldehyde. [Uniprot: P41233]				
Product Application Details					
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin				
Recommended Dilutions	Western Blot 1:1000-1:3000, Immunohistochemistry 1:500, Immunohistochemistry-Paraffin 1:500. Use reported in scientific literature (PMID 19718435)				
Application Notes	In Western blot, a band is observed at ~250 kDa. ABCA1 is known to aggregate. Do not boil samples and use a reducing buffer.				

Images

Western Blot: ABCA1 Antibody (3A1.891.3) [NB400-164] - Western blot analysis of lipid transporters in macrophages treated with or without oxLDL using antibodies against ABCA1 (NB400-164), SR-B1 (NB400-104), and CD36 (NB400-144) . Lrp1Y63F impairs Abca1-mediated cholesterol efflux and increases lipid accumulation in macrophages. Image collected and cropped by CiteAb from the following publication (https://elifesciences.org/articles/29292) licensed under a CC-BY license.

			Ld	Ir+			L	rp1	Y63	; L	dlr"	1-
oxLDL	-	-	-	+	+	+	۲	+	-	+	+	+
Lane	1	2	3	4	5	6	7	8	9	10	11	12
Abca1	-	**	i.	-	-	-	-	-	-		-	
Abcg1	-	-		-	-	-	-	-	-	-	-	-
Cd36	***	*10			-	h	-	**	-	-	-	-
Cub.4	-	-	-	-	-	-	-	-	-	-	-	-
SIDT												



Immunohistochemistry-Paraffin: ABCA1 Antibody (3A1.891.3) [NB400- 164] - Analysis of FFPE mouse liver using ABCA1 antibody at 1:500 on a Bond Rx autostainer (Leica Biosystems). The assay involved 20 minutes of heat induced antigen retrieval (HIER) using 10mM sodium citrate buffer (pH 6.0) and endogenous peroxidase quenching with peroxide block. The sections were incubated with primary antibody for 30 minutes and Bond Polymer Refine Detection (Leica Biosystems) with DAB was used for signal development followed by counterstaining with hematoxylin. Whole slide scanning and capturing of representative images was performed using Aperio AT2 (Leica Biosystems). Plasma membrane staining was observed. Staining was performed by Histowiz.	
Western Blot: ABCA1 Antibody (3A1.891.3) [NB400-164] - Detection of ABCA1 using NB400-164 at various dilutions in ABCA1 transfected HeLa lysates.	1:100 1:200 1:1,600 1:1,600 1:1,5,400 1:1,2,800
Western Blot: ABCA1 Antibody (3A1.891.3) [NB400-164] - Analysis of ABCA1 in (A) T09 treated RAW 264.7 cells and (B) untreated RAW 264.7 cells.	A B 460 268 238 171 117 71 55
Western Blot: ABCA1 Antibody (3A1.891.3) [NB400-164] - Lrp1Y63F mutation downregulates hepatic Abca1 levels.(A) Western blot analysis of Abca1 expression in the liver from Lrp1Y63F;Ldlr-/- & Ldlr-/- mice fed with HCHF diet for 16 weeks (n = 5, each group). (B) Western blot analysis of Abca1 expression in the liver from 16 week HCHF-fed Lrp1Y63F;Ldlr-/- & Ldlr-/- mice after BMT (n = 5, each group). All data are mean ±SEM. **p<0.01. Image collected & cropped by CiteAb from	A $Abca1$





Page 3 of 5 v.20.1 Updated 10/23/2024





Page 4 of 5 v.20.1 Updated 10/23/2024

С

Western Blot: ABCA1 Antibody (3A1.891.3) [NB400-164] - Lrp1Y63F mutation has no effect on atherogenic signaling pathways & morphological phenotype in SMCs.(A) Left: representative images of H&E stained descending thoracic aortas (scale bar = $50 \mu m$). Right: aortic wall thickness & cell number were quantified for the indicated genotypes. Data are expressed mean ±SEM; n = 3 per group & genotype. (B) HCHF diet in Lrp1Y63F;Ldlr-/- mice has no effect on PDGF & TGFβ-mediated signaling in the aorta. Proteins were extracted from whole aorta in Lrp1Y63F;Ldlr-/- & Ldlr-/- mice fed with indicated diets. Expression of proteins involved in PDGF signaling pathways was determined by Western blot analysis. (C) Lrp1Y63F mutation did not affect atherogenic signaling pathways mediated by PDGF in vitro. Primary SMCs from WT & Lrp1Y63F mice were incubated with or without 10 ng/ml PDGF-BB for 10 min. Expression of proteins involved in PDGF signaling were determined by Western blot analysis. (D) Lrp1Y63F mutation did not alter SMC functions in vitro. Primary cultured SMCs from WT & Lrp1Y63F mice were incubated with or without 10 ng/ml PDGF-BB under indicated conditions, then cell migration (left) & proliferation (right) were determined. Values from three independent experiments are expressed as mean ±SEM. *p<0.05. Image collected & cropped by CiteAb from the following publication (https://elifesciences.org/articles/29292), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Ishimaru K, Yoshioka K, Kano K, et al. Sphingosine kinase-2 prevents macrophage cholesterol accumulation and atherosclerosis by stimulating autophagic lipid degradation Sci Rep 2019-12-04 [PMID: 31797978] (WB, Mouse)

Quiroz A, Molina P, Santander N et al. Ovarian Cholesterol Efflux: ABC Transporters and Follicular Fluid HDL Regulate Cholesterol Content in Mouse Oocytes Biol. Reprod. 2019-08-19 [PMID: 31423535]

Xian X, Ding Y, Dieckmann M et al. LRP1 integrates murine macrophage cholesterol homeostasis and inflammatory responses in atherosclerosis Elife. 2017-11-16 [PMID: 29144234] (WB, Mouse)

Sultana A, Cochran BJ, Tabet F et al. Inhibition of inflammatory signaling pathways in 3T3-L1 adipocytes by apolipoprotein A-I. FASEB J. 2016-03-10 [PMID: 26965683] (WB, Mouse)

Do TM, Dodacki A, Alata W et al. Age-Dependent Regulation of the Blood-Brain Barrier Influx/Efflux Equilibrium of Amyloid-beta Peptide in a Mouse Model of Alzheimer's Disease (3xTg-AD). J Alzheimers Dis 2015-01-01 [PMID: 26484906] (Mouse)

Zhou L, Choi HY, Li WP et al. LRP1 controls cPLA2 phosphorylation, ABCA1 expression and cellular cholesterol export. PLoS One. 2009-08-31 [PMID: 19718435] (WB, IHC-P, Mouse)

Do TM, Noel-Hudson MS, Ribes S et al. ABCG2- and ABCG4-mediated efflux of amyloid-beta peptide 1-40 at the mouse blood-brain barrier J Alzheimers Dis 2012-01-01 [PMID: 22391220]

Wang F, Okamoto Y, Inoki I et al. Sphingosine-1-phosphate receptor-2 deficiency leads to inhibition of macrophage proinflammatory activities and atherosclerosis in apoE-deficient mice. J Clin Invest. 120(11):3979-95. doi: 10.1172/JCI42315. 2010-11-01 [PMID: 20978351]

Do TM, Ouellet M, Calon F et al. Direct evidence of abca1-mediated efflux of cholesterol at the mouse blood-brain barrier. Mol Cell Biochem. 2011-06-10 [PMID: 21660464] (WB, Mouse)





Novus Biologicals USA

10730 E. Briarwood Avenue Centennial, CO 80112 USA Phone: 303.730.1950 Toll Free: 1.888.506.6887 Fax: 303.730.1966 nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave Toronto, ON M8Z 4E6 Canada Phone: 905.827.6400 Toll Free: 855.668.8722 Fax: 905.827.6402 canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane Abingdon Science Park Abingdon, OX14 3NB, United Kingdom Phone: (44) (0) 1235 529449 Free Phone: 0800 37 34 15 Fax: (44) (0) 1235 533420 info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com Technical Support: nb-technical@biotechne.com Orders: nb-customerservice@bio-techne.com General: novus@novusbio.com

Products Related to NB400-164

F0105B	Goat anti-Rat IgG Secondary Antibody [Phycoerythrin]
HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
NBP1-30159	Raw 264.7 Whole Cell - T0901317 treated Assay Kit

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB400-164

Earn gift cards/discounts by submitting a publication using this product: www.novusbio.com/publications

www.novusbio.com

