

# Product Datasheet

## Survivin [p Thr34] Antibody [HRP] NB500-236H

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB500-236H](http://www.novusbio.com/NB500-236H)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB500-236H](http://www.novusbio.com/reviews/destination/NB500-236H)



**NB500-236H**

Survivin [p Thr34] Antibody [HRP]

| Product Information         |                                                                                                                                                                               |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                        |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                   |
| Storage                     | Store at 4C in the dark.                                                                                                                                                      |
| Clonality                   | Polyclonal                                                                                                                                                                    |
| Preservative                | No Preservative                                                                                                                                                               |
| Isotype                     | IgG                                                                                                                                                                           |
| Conjugate                   | HRP                                                                                                                                                                           |
| Purity                      | Immunogen affinity purified                                                                                                                                                   |
| Buffer                      | PBS                                                                                                                                                                           |
| Target Molecular Weight     | 16 kDa                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                               |
| Host                        | Rabbit                                                                                                                                                                        |
| Gene ID                     | 332                                                                                                                                                                           |
| Gene Symbol                 | BIRC5                                                                                                                                                                         |
| Species                     | Human, Mouse                                                                                                                                                                  |
| Reactivity Notes            | Human and mouse reactivity reported in scientific literature (PMID: 17510430).                                                                                                |
| Specificity/Sensitivity     | Survivin [p Thr34] Antibody [NB500-236] is specific for phosphorylated survivin.                                                                                              |
| Immunogen                   | This Survivin [p Thr34] Antibody was developed against a synthetic peptide with a phosphorylated Threonine (amino acid 34) corresponding to human Survivin. [UniProt# O15392] |
| Product Application Details |                                                                                                                                                                               |
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence                                                                                                                         |
| Recommended Dilutions       | Western Blot, Immunocytochemistry/ Immunofluorescence                                                                                                                         |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                        |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NB500-236H**

---

|                |                                           |
|----------------|-------------------------------------------|
| NBP2-24891H    | Rabbit IgG Isotype Control [HRP]          |
| NB500-236G     | Survivin [p Thr34] Antibody [DyLight 488] |
| NB110-92717PEP | Survivin Antibody Blocking Peptide        |
| DVE00          | VEGF [HRP]                                |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB500-236H](http://www.novusbio.com/reviews/submit/NB500-236H)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

