

# Product Datasheet

## Survivin Antibody (60.11) [DyLight 650] NB500-238C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB500-238C](http://www.novusbio.com/NB500-238C)

Updated 11/11/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB500-238C](http://www.novusbio.com/reviews/destination/NB500-238C)



**NB500-238C**

Survivin Antibody (60.11) [DyLight 650]

| Product Information     |                                                                                             |
|-------------------------|---------------------------------------------------------------------------------------------|
| Unit Size               | 0.1 ml                                                                                      |
| Concentration           | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage                 | Store at 4C in the dark.                                                                    |
| Clonality               | Monoclonal                                                                                  |
| Clone                   | 60.11                                                                                       |
| Preservative            | 0.05% Sodium Azide                                                                          |
| Isotype                 | IgG2a Kappa                                                                                 |
| Conjugate               | DyLight 650                                                                                 |
| Purity                  | Protein G purified                                                                          |
| Buffer                  | 50mM Sodium Borate                                                                          |
| Target Molecular Weight | 16 kDa                                                                                      |

| Product Description     |                                                                                                                             |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Host                    | Mouse                                                                                                                       |
| Gene ID                 | 332                                                                                                                         |
| Gene Symbol             | BIRC5                                                                                                                       |
| Species                 | Human, Mouse, Rat                                                                                                           |
| Specificity/Sensitivity | Survivin Antibody (60.11) [NB500-238] is specific for the cytoplasmic form of Survivin.                                     |
| Immunogen               | This Survivin Antibody (60.11) [DyLight 650] was developed against full length recombinant human Survivin [UniProt# O15392] |
| Notes                   | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.                                           |

| Product Application Details |                                                                                                                                                                                                                                     |
|-----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry-Paraffin, ELISA, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Proximity Ligation Assay, Knockdown Validated |
| Recommended Dilutions       | Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Proximity Ligation Assay, Flow (Intracellular), Knockdown Validated |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                              |

**Publications**

Kusner LL, Ciesielski MJ, Marx A et al. Survivin as a potential mediator to support autoreactive cell survival in myasthenia gravis: a human and animal model study. PLoS ONE 2014-07-23 [PMID: 25050620] (FLOW, Human)

Details:

Using the DyLight 650 conjugated version of NB500-238, catalog number NB500-238C.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB500-238C**

|                |                                                        |
|----------------|--------------------------------------------------------|
| NBP1-96981C    | Mouse IgG2a Kappa Isotype Control (M2AK) [DyLight 650] |
| NB500-238F     | Survivin Antibody (60.11) [FITC]                       |
| NB110-92717PEP | Survivin Antibody Blocking Peptide                     |
| DVE00          | VEGF [HRP]                                             |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB500-238C](http://www.novusbio.com/reviews/submit/NB500-238C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

